

# Product Datasheet

## Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free NBP2-80614

Unit Size: 0.2 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-80614](http://www.novusbio.com/NBP2-80614)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-80614](http://www.novusbio.com/reviews/destination/NBP2-80614)




**NBP2-80614**

Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free


Product Information	
Unit Size	0.2 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CPP32 4-1-18
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	31.7 kDa
Product Description	
Description	Novus Biologicals Mouse Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free (NB500-210) is a monoclonal antibody validated for use in IHC, WB, ICC/IF, Simple Western and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	836
Gene Symbol	CASP3
Species	Human, Mouse, Rat
Immunogen	This Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free was developed against full-length recombinant human Caspase 3 [UniProt# P42574].
Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:1000, Simple Western 1:2000, Immunohistochemistry 1:200-1:500, Immunocytochemistry/ Immunofluorescence 2 ug/ml, Immunoprecipitation 2 ug / mg lysate, Immunohistochemistry-Paraffin 1:200-1:500, Immunohistochemistry-Frozen reported in scientific literature (PMID 31242448)
Application Notes	<p>A band is seen at ~32 kDa for the inactive form of Caspase 3 and ~17-22 kDa for the active form of Caspase 3 in Western Blot. (See protocol for additional information).</p> <p>In Simple Western only 10 - 15 uL of the recommended dilution is used per data point.</p> <p>See <a href="#">Simple Western Antibody Database</a> for Simple Western validation: antibody dilution of 1:2000. Separated by Size-Wes, Sally Sue/Peggy Sue.</p>

## Images

Western Blot: Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free [NBP2-80614] - Detection of Caspase (19 and 35 kDa) from HEK293 cell extract using Caspase-3 Antibody (CPP32 4-1-18) (NB500-210). Lanes 1 and 2 contain inactive and active Caspase, respectively. Image from the standard format of this antibody.




Immunocytochemistry: Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free [NBP2-80614] - Caspase-3 was detected in immersion fixed paraffin-embedded sections of human bladder 1:300 dilution of mouse monoclonal Caspase-3 Antibody (CPP32 4-1-18) (NB500-210, Novus Biologicals), for 1 hour at room temperature followed by anti-mouse IgG VisUCyte H




Immunohistochemistry-Paraffin: Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free [NBP2-80614] - Rat epithelial cells of the tongue base. Antigen retrieval method: Citrate buffer. Image from the standard format of this antibody.



Immunohistochemistry: Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free [NBP2-80614] - Analysis of a FFPE human spleen section using 1:200 dilution of Caspase-3 Antibody (CPP32 4-1-18). The staining was developed using HRP conjugated anti-mouse secondary antibody and DAB reagent. This Caspase 3 antibody generated a specific staining in the



Simple Western: Caspase-3 Antibody (CPP32 4-1-18) - Azide and BSA Free [NBP2-80614] - Lane view shows a specific band for Caspase 3 in 0.5 mg/ml of Hek293 lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system. Image from the standard format of this antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-80614**

---

NBP3-11853	Jurkat Staurosporine Treated / Untreated Cell Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-80614](http://www.novusbio.com/reviews/submit/NBP2-80614)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

