

Product Datasheet

EGLN1/PHD2 Antibody (2445D) - Azide and BSA Free NBP2-80714

Unit Size: 0.1 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80714

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80714



NBP2-80714

EGLN1/PHD2 Antibody (2445D) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	2445D
Preservative	No Preservative
Purity	Protein A or G purified
Buffer	PBS
Target Molecular Weight	46 kDa


Product Description	
Description	Novus Biologicals Rabbit EGLN1/PHD2 Antibody (2445D) - Azide and BSA Free (NBP2-78782) is a recombinant monoclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	54583
Gene Symbol	EGLN1
Species	Human, Mouse
Immunogen	This EGLN1/PHD2 antibody was developed against <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human EGLN1/PHD2 protein (amino acids 2-426) [UniProt #Q9QZT9].

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/mL, Immunohistochemistry 3 - 25 ug/mL, Immunohistochemistry-Paraffin 3 - 25 ug/mL




Images

Western Blot: EGLN1/PHD2 Antibody (2445D) - Azide and BSA Free [NBP2-80714] - Western blot shows lysates of PC3 human prostate cancer cell line, HEK293T human embryonic kidney cell line, and NIH3T3 mouse embryonic fibroblast cell line. PVDF membrane was probed with 1 ug/mL of Rabbit Anti-Human/Mouse EGLN1/PHD2 Monoclonal Antibody (Catalog # NBP2-78782) followed by HRP-conjugated Anti-Rabbit IgG Secondary Antibody (Catalog # HAF008). A specific band was detected for EGLN1/PHD2 at approximately 50 kDa (as indicated). This experiment was conducted under reducing conditions. Image from the standard format of this antibody.



Immunohistochemistry: EGLN1/PHD2 Antibody (2445D) - Azide and BSA Free [NBP2-80714] - EGLN1/PHD2 was detected in immersion fixed paraffin-embedded sections of human kidney using Rabbit Anti-Human/Mouse EGLN1/PHD2 Monoclonal Antibody (Catalog # NBP2-78782) at 3 ug/mL for 1 hour at room temperature followed by incubation with an anti-rabbit secondary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80714

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
H00054583-P01-10ug	Recombinant Human EGLN1/PHD2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80714

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

