

Product Datasheet

FGFR1 Antibody (M2F12) - Azide and BSA Free NBP2-80734

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80734

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80734



NBP2-80734


FGFR1 Antibody (M2F12) - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	M2F12
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein G purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Mouse FGFR1 Antibody (M2F12) - Azide and BSA Free (NB120-829) is a monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2260
Gene Symbol	FGFR1
Species	Human, Mouse, Rat
Immunogen	Recombinant human ectodomain of FGFR1a expressed in E. coli beginning with pro23; antigen contained NH2-terminal gly-ser-pro-gly-ile and COOH-terminal glu-phe sequences.
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, CyTOF-ready
Recommended Dilutions	Western Blot 1 ug/ml, Flow Cytometry 1:500 - 1:1000, Immunohistochemistry reported in scientific literature (PMID 8168112), Immunocytochemistry/ Immunofluorescence 10 ug/ml, Immunoprecipitation 1:10 - 1:500, Immunohistochemistry-Frozen 1:10 - 1:500, CyTOF-ready
Application Notes	Western blotting: 1 ug/mL recognizes 10 ng recombinant bacterial, baculo-viral, or native FGFR1. Immunoprecipitation: Native, nuclear, recombinant, and ligand-labeled FGFR1. Customers have reported difficulties using this antibody in immunocytochemistry staining (see review). This antibody does not appear to work in IHC paraffin embedded sections.




Images


Western Blot: FGFR1 Antibody (M2F12) - Azide and BSA Free [NBP2-80734] - Total protein from Ntera-2 was separated on a 7.5% gel by SDS-PAGE. Protein was transferred to PVDF membrane, blocked and probed with 1 ug/ml of anti-FGFR1. FGFR1 protein was detected at approximately 150 kDa using an anti-mouse HRP secondary antibody. Image from the standard format of this antibody.




Immunocytochemistry/Immunofluorescence: FGFR1 Antibody (M2F12) - Azide and BSA Free [NBP2-80734] - NTERA2 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton-X100. The cells were incubated with FGFR1 (M2F12) (NB120-829) at a 1:100 dilution overnight at 4 degrees Celsius and detected with Dylight 488 (Green) at a 1:500 dilution. Actin was detected with Phalloidin 568 (Red) at a 1:200 dilution. Nuclei were detected with DAPI (Blue) at 2.0 ug/ml in 1X PBS. Cells were imaged using a 40X objective. Image from the standard format of this antibody.



Flow Cytometry: FGFR1 Antibody (M2F12) - Azide and BSA Free [NBP2-80734] - Analysis of FGF R1 in HEK293 cells (2×10^6) stained with FGF R1 antibody (NB120-829, red) at a 1:1000 dilution and detected with a FITC conjugated goat anti-mouse IgG. Shown with mIgG2a isotype control (green). Two distinct samples shown. Image from the standard format of this antibody.



Flow Cytometry: FGFR1 Antibody (M2F12) - Azide and BSA Free [NBP2-80734] - Analysis of FGF R1 in MCF-7 cells (1×10^6 cells/ml) stained with FGF r1 antibody (NB120-829, red peaks) at a 1:1000 dilution. Detected with a FITC conjugated goat anti-mouse IgG. Shown with mIgG2a isotype control (green). Two distinct samples shown. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80734

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80734

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

