

Product Datasheet

IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free NBP2-80807

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80807

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80807



NBP2-80807

IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	34N4G11
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Mouse IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free (NBP2-19015) is a monoclonal antibody validated for use in IHC and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9466
Gene Symbol	IL27RA
Species	Human
Immunogen	A portion of amino acids 200-250 of human IL-27R (WSX-1) was used as the immunogen for this antibody.


Product Application Details	
Applications	Immunohistochemistry-Paraffin, Flow Cytometry, Flow (Cell Surface), Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Flow Cytometry 1 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin 5 ug/ml, Flow (Cell Surface) 1 ug/ml, CyTOF-ready
Application Notes	This antibody is Cytof ready.

Images


Immunohistochemistry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - IHC-P detection of IL27RA protein in a section of normal human colon tissue using IL27RA antibody (34N4G11) at a concentration of 5ug/ml. Intense membrane with cytoplasmic staining was observed in the cells of crypts and the absorptive columnar epithelium




Flow (Cell Surface): IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - Analysis using the Alexa Fluor (R) 647 conjugate of NBP2-19015. Surface staining of IL27R on freshly isolated (PPI negative) human PBMCs using IL-27R antibody and CD3 antibody at 10 ul/10⁶ cells. Flow kit, mouse IgG2b isotype control antibody, and RBC lysis buffer were used for this test.




Immunohistochemistry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - Adenocarcinoma of the stomach probed with IL-27R antibody at 5 ug/ml. Image from the standard format of this antibody.




Immunohistochemistry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - FFPE human prostate probed with IL-27R antibody at 5 ug/ml. Image from the standard format of this antibody.




Immunohistochemistry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - IHC-P detection of IL27RA protein in a section of human colo-rectal cancer tissue using IL27RA antibody (34N4G11) at a concentration of 5ug/ml. The cancer cells showed specific positivity for IL27RA and the staining intensity was more towards the cellular




Immunohistochemistry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - IHC-P detection of IL27RA protein in a section of human renal/kidney cancer tissue using IL27RA antibody (34N4G11) at a concentration of 5ug/ml. The representative image shows membrane-cytoplasmic positivity for IL27RA protein in the cancer cells. Image f



Immunohistochemistry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - IHC-P detection of IL27RA protein in a section of human stomach cancer tissue using IL27RA antibody (34N4G11) at a concentration of 5ug/ml. The cancer cells showed intense membrane-cytoplasmic positivity for IL27RA protein, whereas, the glands in adjacent



Flow Cytometry: IL-27R alpha/WSX-1/TCCR Antibody (34N4G11) - Azide and BSA Free [NBP2-80807] - Surface staining on human peripheral blood monocytes using IL-27R antibody at 1 ug/10⁶ cells (red represents specific antibody and green represents isotype control). APC-conjugated goat anti-mouse Ig F(ab')₂ secondary antibody, CS flow kit and an FcR block were used for this test. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80807

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP3-21361PEP	IL-27R alpha/WSX-1/TCCR Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80807

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

