

Product Datasheet

MGMT Antibody (MT 23.2) - Azide and BSA Free NBP2-80847

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80847

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80847



NBP2-80847

MGMT Antibody (MT 23.2) - Azide and BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	MT 23.2
Preservative	No Preservative
Isotype	IgG2b
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	24 kDa
Product Description	
Description	Novus Biologicals Mouse MGMT Antibody (MT 23.2) - Azide and BSA Free (NB100-168) is a monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and Simple Western. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	4255
Gene Symbol	MGMT
Species	Human, Mouse (Negative)
Reactivity Notes	This antibody is mouse reactivity negative.
Immunogen	Recombinant human MGMT (O6-methylguanine-DNA methyltransferase) purified from E. coli [UniProt# P16455]
Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, CyTOF-ready
Recommended Dilutions	Western Blot 1:1000, Simple Western 1:1000, Flow Cytometry 1-2 ug/ml, Immunohistochemistry 1:500-1000, Immunocytochemistry/ Immunofluorescence 1:500, Immunohistochemistry-Paraffin 1:500-1:1000, Flow (Intracellular), CyTOF-ready
Application Notes	<p>This MGMT (MT 23.2) antibody is useful for Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry on paraffin-embedded sections and Western Blot. In WB, a band at approx. 24 kDa can be seen. In ICC/IF, strong staining has been seen in CEM cells and no signal has been seen in TK6 cells.</p> <p>In Simple Western only 10 - 15 ul of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: antibody dilution of 1:1000. Separated by Size-Wes, Sally Sue/Peggy Sue.</p>

Images


Western Blot: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - Detection of MGMT in MCF7 lysate using NB100-168 at 1:1000 dilution. Image from the standard format of this antibody.

250>
150>
100>
75>
50>
37>
25>
20>
15>
10>


Immunocytochemistry/Immunofluorescence: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - MCF7 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton-X100. The cells were incubated with anti-MGMT (MT 23.2) NB100-168 at a 1:200 dilution overnight at 4C and detected with an anti-mouse Dylight 488 (Green) at a 1:500 dilution. Actin was counterstained with Phalloidin 568 (Red) at a 1:200 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective. Image from the standard format of this antibody.




Immunohistochemistry-Paraffin: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - IHC analysis of a formalin fixed paraffin-embedded (FFPE) human kidney using 1:2000 conc. of MGMT antibody (clone MT 23.2) on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of heat induced antigen retrieval (HIER) using 10mM sodium




Flow Cytometry: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - Detection of MGMT in Human Jurkat Cell Line by Flow Cytometry. Human Jurkat cell line was stained with Mouse Anti- MGMT Monoclonal Antibody (Catalog # NB100-168, filled histogram), or Mouse IgG2B isotype control (Catalog # MAB0041, open histogram) followed by APC-conjugated Anti-Mouse IgG Secondary Antibody (Catalog # F0101B). To facilitate intracellular staining, cells were fixed with FlowX FoxP3 Fixation & Permeabilization Buffer Kit (Catalog # FC012). Images may not be copied, printed or otherwise disseminated without express written permission of Novus Biologicals a Bio-Techne brand. Image from the standard format of this antibody.



Immunocytochemistry/Immunofluorescence: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - Detection of MGMT (Green) in HeLa cells using NB100-168 at a 1:50 dilution. Nuclei (Blue) are counterstained using Hoechst 33258. Image from the standard format of this antibody.




Flow Cytometry: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - An intracellular stain was performed on HeLa cells with MGMT Antibody (MT 23.2) antibody NB100-168 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 1 ug/mL for 30 minutes at room temperature, followed by mouse F(ab)2 IgG (H+L) PE-conjugated secondary antibody (F0102B, R&D Systems). Image from the standard format of this antibody.



Copyright © 2018 Novus Biologicals

Simple Western: MGMT Antibody (MT 23.2) - Azide and BSA Free [NBP2-80847] - Simple Western lane view shows a specific band for MGMT in 0.2 mg/ml of MCF-7 lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80847

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80847

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

