

Product Datasheet

Transportin 1 Antibody (D45) - Azide and BSA Free NBP2-80991

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80991

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80991




NBP2-80991**Transportin 1 Antibody (D45) - Azide and BSA Free**

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	D45
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Mouse Transportin 1 Antibody (D45) - Azide and BSA Free (NB600-1397) is a monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3842
Gene Symbol	TNPO1
Species	Human, Mouse, Bovine, Canine, Monkey, Rabbit
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 15738418).
Immunogen	Human recombinant His-tagged Transportin 1. [UniProt# Q92973]
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Microarray, CyTOF-ready
Recommended Dilutions	Western Blot 1:1000, Flow Cytometry, Immunohistochemistry 1:200-1:1000, Immunocytochemistry/ Immunofluorescence reported in scientific literature (PMID 9298975), Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:200-1:1000. Use reported in scientific literature (PMID 24647604), Microarray, Flow (Intracellular), CyTOF-ready
Application Notes	WB: Detects a band of approximately 93 kDa. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors. This antibody is CyTOF ready.




Images


Western Blot: Transportin 1 Antibody (D45) - Azide and BSA Free [NBP2-80991] - Whole cell protein from A431, Hek293 and MCF7 cells was separated on a 7.5% gel by SDS-PAGE, transferred to PVDF membrane and blocked in 5% non-fat milk in TBST. The membrane was probed with 1.0 ug/ml anti-Transportin 1 (D45) in 1% milk, and detected with an anti-mouse HRP secondary antibody using chemiluminescence. Image from the standard format of this antibody.




Immunohistochemistry-Paraffin: Transportin 1 Antibody (D45) - Azide and BSA Free [NBP2-80991] - IHC analysis of a formalin fixed paraffin embedded tissue section of human pancreas using Transportin 1 antibody (clone D45) at 1:200 dilution. The antibody generated specific nuclear and cytoplasmic signal of Transportin 1 protein in the glandular cells



Western Blot: Transportin 1 Antibody (D45) - Azide and BSA Free [NBP2-80991] - Human lung cancer cell line usig lysates from HCC44 cell line. (N: Nuclear fraction, C: Cytoplasmic fraction). Image from verified customer review. Image from the standard format of this antibody.



Immunohistochemistry-Paraffin: Transportin 1 Antibody (D45) - Azide and BSA Free [NBP2-80991] - IHC analysis of a formalin fixed paraffin embedded tissue section of human pancreas using Transportin 1 antibody (clone D45) at 1:200 dilution. The antibody generated specific nuclear and cytoplasmic signal of Transportin 1 protein in the glandular cells





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80991

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80991

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

