

Product Datasheet

IFN-gamma Antibody (D9D10) - Chimeric - Azide and BSA Free NBP2-81044-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81044

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81044



NBP2-81044-0.2mg

IFN-gamma Antibody (D9D10) - Chimeric - Azide and BSA Free

Product Information	
Unit Size	0.2 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	D9D10
Preservative	0.02% Proclin 300
Isotype	IgG Kappa
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit IFN-gamma Antibody (D9D10) - Chimeric - Azide and BSA Free (NBP2-81044) is a recombinant monoclonal antibody validated for use in WB, ELISA and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3458
Gene Symbol	IFNG
Species	Human
Reactivity Notes	Cross Reactivity: Baboon
Specificity/Sensitivity	Binds human IFN-gamma, a dimerized soluble cytokine that is the only member of the type II class of interferons.
Immunogen	FPLC-purified natural human IFN-gamma. Refer to the original article [PMID: 3114127]

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Block/Neutralize
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Block/Neutralize
Application Notes	This chimeric rabbit antibody was made using the variable domain sequences of the original murine IgG1 format, for improved compatibility with existing reagents, assays and techniques. This antibody can be used to detect both native and recombinant intereferon gamma (IFNg). It is reported only to bind to the functional, dimeric form of IFNg.



Images

Western Blot: IFN-gamma Antibody (D9D10) [NBP2-81044] - Western Blot using anti-IFN gamma antibody (D9D10) [NBP2-81044]. Human thymus lysate (35ug protein in RIPA buffer) was resolved on a 15% SDS PAGE gel and blots probed with the rabbit IgG-chimeric version of D9D10 [NBP2-81044] at 0.1 ug/ml before detection using an anti-rabbit secondary antibody. Protein was detected by chemiluminescence. The expected unmodified running size for IFN gamma is 19.3kDa, though the form shown here is likely glycosylated. NBP2-81044 successfully detected IFN gamma in human thymus samples.



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81044-0.2mg

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-34992-100ug	Recombinant Human IFN-gamma Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81044

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

