

Product Datasheet

Mouse Procollagen III C-Terminal ProPeptide ELISA Kit (Colorimetric) NBP2-81203

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81203

Updated 9/9/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81203




NBP2-81203**Mouse Procollagen III C-Terminal ProPeptide ELISA Kit (Colorimetric)**


Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Species	Mouse
Specificity/Sensitivity	This kit recognizes Mouse Procollagen III C-Terminal ProPeptide in samples. No significant cross-reactivity or interference between Mouse Procollagen III C-Terminal ProPeptide and analogues was observed.
Kit Components	Micro ELISA Plate (Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Manual
Standard Curve Range	15.63 - 1000 pg/mL
Sensitivity	9.38 pg/mL
Inter Assay	CV% < 4.21%
Intra Assay	CV% < 5.24%
Assay Type	Sandwich-ELISA
Spiking Recovery	85-107%
Suitable Sample Type	Serum, plasma, cell supernatant, and other biological fluids
Sample Volume	100 ul
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Mouse Procollagen III C-Terminal ProPeptide ELISA Kit (Colorimetric) [NBP2-81203] - Samples were spiked with high concentrations of Mouse Procollagen III C-Terminal ProPeptide and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Mouse Procollagen III C-Terminal ProPeptide ELISA Kit (Colorimetric) [NBP2-81203] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81203

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

