

Product Datasheet

Cyclin T1 Antibody - BSA Free NBP2-81898

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81898

Updated 3/4/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81898



NBP2-81898


Cyclin T1 Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	80 kDa
Product Description	
Host	Rabbit
Gene ID	904
Gene Symbol	CCNT1
Species	Human, Mouse, Rat
Immunogen	Rabbit polyclonal Cyclin T1 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human Cyclin T1. The immunogen is located within the last 50 amino acids of Cyclin T1. Amino Acid Sequence: NTDKKPRPPPLPSEPPPL
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin




Images


Western Blot: Cyclin T1 Antibody [NBP2-81898] - Analysis of Cyclin T1 in rat brain tissue lysate with Cyclin T1 antibody at 1 ug/ml in (A) the absence and (B) presence of blocking peptide.




Immunocytochemistry/ Immunofluorescence: Cyclin T1 Antibody - BSA Free [NBP2-81898] - Immunofluorescence of Cyclin T1 in rat brain tissue with Cyclin T1 antibody at 20 ug/ml.



Immunohistochemistry: Cyclin T1 Antibody [NBP2-81898] - Analysis of Cyclin T1 in rat brain tissue with Cyclin T1 antibody at 2.5 ug/ml.



Immunocytochemistry/ Immunofluorescence: Cyclin T1 Antibody - BSA Free [NBP2-81898] - Immunofluorescence of Cyclin T1 in rat brain tissue with Cyclin T1 antibody at 20 u/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81898

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00000904-Q01-10ug	Recombinant Human Cyclin T1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81898

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

