

Product Datasheet

LMX1A Antibody - BSA Free

NBP2-84141-0.1ml

Unit Size: 0.1ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-84141

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-84141



NBP2-84141-0.1ml

LMX1A Antibody - BSA Free


Product Information	
Unit Size	0.1ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Protein A purified
Buffer	PBS, 2% Sucrose
Product Description	
Description	Novus Biologicals Rabbit LMX1A Antibody - BSA Free (NBP2-84141) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	4009
Gene Symbol	LMX1A
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human LMX1A. Peptide sequence: QNQRAKMKKLARRQQQQQDQQNTQRLSSAQTNGGGSAGMEGIMNPYTAL The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin




Images

Western Blot: LMX1A Antibody [NBP2-84141] - Host: Rabbit. Target Name: LMX1A. Sample Type: MCF7. Antibody Dilution: 1.0ug/ml


LMX1A



Immunohistochemistry: LMX1A Antibody [NBP2-84141] - Rabbit Anti-LMX1A Antibody. Paraffin Embedded Tissue: Human Muscle. Cellular Data: Skeletal muscle cells. Antibody Concentration: 4.0-8.0 ug/ml. Magnification: 400X




Western Blot: LMX1A Antibody [NBP2-84141] - WB Suggested Anti-LMX1A Antibody Titration: 0.3125ug/ml. ELISA Titer: 1:62500. Positive Control: Jurkat cell lysate



Western Blot: LMX1A Antibody [NBP2-84141] - Host: Rabbit. Target Name: TRIM10. Sample Type: 721_B. Antibody Dilution: 1.0ug/ml

LMX1A





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-84141-0.1ml

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
H00004009-P01-10ug	Recombinant Human LMX1A GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-84141

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

