

Product Datasheet

hnRNP AB Antibody - BSA Free NBP2-87582

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-87582

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-87582



NBP2-87582

hnRNP AB Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Protein A purified
Buffer	PBS, 2% Sucrose

Product Description	
Description	Novus Biologicals Rabbit hnRNP AB Antibody - BSA Free (NBP2-87582) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3182
Gene Symbol	HNRNPAB
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of human hnRNP AB. Peptide sequence: QQGYGPGYGGYDYSPTYGYYGYGPGYDYSQGSTNYGKSQRRGGHQNNYKPY The peptide sequence for this immunogen was taken from within the described region.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin



Images


Western Blot: hnRNP AB Antibody [NBP2-87582] - WB Suggested Anti-HNRPAB Antibody Titration: 1.25ug/ml. Positive Control: Daudi cell lysate. HNRPAB is supported by BioGPS gene expression data to be expressed in Daudi




Immunohistochemistry-Paraffin: hnRNP AB Antibody [NBP2-87582] - Rabbit Anti-HNRPAB Antibody. Paraffin Embedded Tissue: Human cardiac cell. Cellular Data: Epithelial cells of renal tubule. Antibody Concentration: 4.0-8.0 ug/ml. Magnification: 400X



Immunohistochemistry-Paraffin: hnRNP AB Antibody [NBP2-87582] - Rabbit Anti-HNRPAB Antibody. Paraffin Embedded Tissue: Human Lung. Cellular Data: Alveolar cells. Antibody Concentration: 4.0-8.0 ug/ml. Magnification: 400X



Immunohistochemistry-Paraffin: hnRNP AB Antibody [NBP2-87582] - Rabbit Anti-HNRPAB Antibody. Paraffin Embedded Tissue: Human Heart. Cellular Data: Myocardial cells. Antibody Concentration: 4.0-8.0 ug/ml. Magnification: 400X





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-87582

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-23119	Recombinant Human hnRNP AB His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-87582

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

