

Product Datasheet

DLX3 Antibody - BSA Free NBP2-88803

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-88803

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-88803



NBP2-88803

DLX3 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Protein A purified
Buffer	PBS, 2% Sucrose


Product Description	
Description	Novus Biologicals Rabbit DLX3 Antibody - BSA Free (NBP2-88803) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1747
Gene Symbol	DLX3
Species	Human, Mouse
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of human DLX3. Peptide sequence: SASPSYLDDPTNSWYHAQNLSGPHLQQQPPQPATLHHASPGPPPNGAVY The peptide sequence for this immunogen was taken from within the described region.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml



Images


Western Blot: DLX3 Antibody [NBP2-88803] - Host: Mouse. Target Name: DLX3. Sample Tissue: Mouse Liver. Antibody Dilution: 1ug/ml



Western Blot: DLX3 Antibody [NBP2-88803] - WB Suggested Anti-DLX3 Antibody Titration: 1.25ug/ml. ELISA Titer: 1:1562500. Positive Control: Transfected 293T



Western Blot: DLX3 Antibody [NBP2-88803] - Host: Rabbit. Target: DLX3. Positive control (+): MCF7 (N10). Negative control (-): HepG2 (HG). Antibody concentration: 1ug/ml





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-88803

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-23004	Recombinant Human DLX3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-88803

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

