

Product Datasheet

Human Coronavirus Spike Protein Antibody (028) - Azide and BSA Free NBP2-90334

Unit Size: 50 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-90334

Updated 2/26/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-90334



NBP2-90334**Human Coronavirus Spike Protein Antibody (028) - Azide and BSA Free**

Product Information	
Unit Size	50 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	028
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	0.2 um filtered solution in PBS
Product Description	
Description	This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Gene ID	3200426
Gene Symbol	S
Species	Virus
Reactivity Notes	Human coronavirus (HCoV-HKU1)
Specificity/Sensitivity	HCoV-HKU1 (isolate N1) Spike/S1 No cross-reactivity in ELISA with: HCoV-HKU1 (isolate N5) Spike S1+S2 ECD Protein SARS-CoV-2 Spike S1 Protein SARS-CoV-2 Spike S1 NTD Protein SARS-CoV-2 Spike S1 NTD Protein SARS-CoV-2 Spike RBD Protein SARS-CoV-2 Spike S2 Protein SARS-CoV-2 Spike S1+S2 ECD Protein SARS-CoV Spike S1+S2 ECD (S577A, Isolate Tor2) Protein MERS-CoV Spike S1+S2 ECD Protein HCoV-NL63 Spike S1+S2 ECD Protein HCoV-229E Spike S1+S2 ECD Protein HCoV-OC43 Spike S1+S2 ECD Protein
Immunogen	This antibody was obtained from a rabbit immunized with purified, recombinant HCoV-HKU1 (Accession#: YP_173238.1; Met1-Arg760).
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA 1:5000-1:10000
Application Notes	Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-90334

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-90334

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

