

# Product Datasheet

## S100A1 Antibody (049) [PerCP] NBP2-90504PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-90504PCP](http://www.novusbio.com/NBP2-90504PCP)

Updated 7/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-90504PCP](http://www.novusbio.com/reviews/destination/NBP2-90504PCP)



**NBP2-90504PCP**

S100A1 Antibody (049) [PerCP]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 049   |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG   |
| Conjugate           | PerCP   |
| Purity              | Protein A purified  |
| Buffer              | PBS   |

| Product Description |  |
|---------------------|--|
| Host                | Rabbit   |
| Gene ID             | 6271   |
| Gene Symbol         | S100A1   |
| Species             | Mouse  |
| Immunogen           | This antibody was obtained from a rabbit immunized with purified, recombinant Mouse S100A1 (Uniprot#: P56565; Gly2-Ser94). |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | ELISA, Sandwich ELISA Capture  |
| Recommended Dilutions       | ELISA, Sandwich ELISA Capture  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-90504PCP**

---

|                  |                                      |
|------------------|--------------------------------------|
| NBP2-24891PCP    | Rabbit IgG Isotype Control [PerCP]   |
| NBP1-72335-0.1mg | Recombinant Mouse S100A1 His Protein |
| 9705-S1-100      | S100A1 [Unconjugated]                |
| AF1356           | CD117/c-kit Antibody [Unconjugated]  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-90504PCP](http://www.novusbio.com/reviews/submit/NBP2-90504PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

