

Product Datasheet

CYLD Antibody - BSA Free

NBP2-92436-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-92436

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-92436



NBP2-92436-0.1ml

CYLD Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	107 kDa

Product Description	
Description	Novus Biologicals Rabbit CYLD Antibody - BSA Free (NBP2-92436) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1540
Gene Symbol	CYLD
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 654-953 of human CYLD (NP_001035877.1). TKIMKLRKILEKVEAASGFTSEEKDP EEFLNLFHHILRVEPLLKIRSAGQKVQDC YFYQIFMEKNEKVGVP TIQQLLEWSFINSLKFAEAPSLIIQM PRFGKDFKLFK KIFPSLELNITDLLEDTPRQCRICGGLAMYECRECYDDPDISAGKIKQFCKTCNT QVHLHPKRLN HKYNPVSLPKDLPDWDWRHGCIPCQNMELFAVLCIETSHYVAF VKYGGKDDSAWLFFDSMADRGGQNGFNIPQVTPCPEVGEYLKMSLEDLHSLD SRRIQGCARRLLCDA YMCMYQSPTMSLYK


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:100 - 1:500, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Paraffin

Images


Western Blot: CYLD Antibody [NBP2-92436] - Western blot analysis of extracts of various cell lines, using CYLD antibody (NBP2-92436) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.




Immunohistochemistry-Paraffin: CYLD Antibody [NBP2-92436] - Immunohistochemistry of paraffin-embedded mouse spleen using CYLD Rabbit pAb (NBP2-92436) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: CYLD Antibody [NBP2-92436] - Immunohistochemistry of paraffin-embedded rat brain using CYLD Rabbit pAb (NBP2-92436) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.




Immunohistochemistry-Paraffin: CYLD Antibody [NBP2-92436] - Immunohistochemistry of paraffin-embedded rat spleen using CYLD Rabbit pAb (NBP2-92436) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry-Paraffin: CYLD Antibody [NBP2-92436] - Immunohistochemistry of paraffin-embedded human lung cancer using CYLD Rabbit pAb (NBP2-92436) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry-Paraffin: CYLD Antibody [NBP2-92436] -
Immunohistochemistry of paraffin-embedded human esophageal cancer
using CYLD Rabbit pAb (NBP2-92436) at dilution of 1:100 (40x lens).
Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0
before commencing with IHC staining protocol.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-92436-0.1ml

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-92436

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

