

Product Datasheet

Protein Kinase A regulatory subunit I alpha Antibody - BSA Free NBP2-93193-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-93193

Updated 2/17/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-93193



NBP2-93193-0.1ml

Protein Kinase A regulatory subunit I alpha Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol

Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit Protein Kinase A regulatory subunit I alpha Antibody - BSA Free (NBP2-93193) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5573
Gene Symbol	PRKAR1A
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Protein Kinase A regulatory subunit I alpha (NP_002725.1). MESGTAASEEARSLRECELYVQKHNIQALLKDSIVQLCTARPERPMAFLREYF ERLEKEEAKQIQNLQKAGTRTDSREDEISPPPPNPVVKGRRRRGAISAEVYTEE DAASYVRKVIPKDYKTMAALAKAIEKNVLFSLDDNERSDIFDAMFSVSFIAGET VIQQGDEGDNFYVIDQGETDVYVNNNEWATSVGEGGSFGELALIYGTPRAATVK AKTNVKLWGIDRDSYRRILMGSTLRKRKMYEEFL


Product Application Details	
Applications	Western Blot, Knockout Validated
Recommended Dilutions	Western Blot 1:500 - 1:1000, Knockout Validated

Images


Western Blot: Protein Kinase A regulatory subunit I alpha Antibody [NBP2-93193] - Western blot analysis of extracts from normal (control) and Protein Kinase A regulatory subunit I alpha knockout (KO) 293T cells, using Protein Kinase A regulatory subunit I alpha antibody (NBP2-93193) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 180s.




Western Blot: Protein Kinase A regulatory subunit I alpha Antibody [NBP2-93193] - Analysis of extracts of various cell lines, using Protein Kinase A regulatory subunit I alpha at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.




Western Blot: Protein Kinase A regulatory subunit I alpha Antibody - Azide and BSA Free [NBP2-93193] - Western blot analysis of lysates from Mouse testis using Protein Kinase A regulatory subunit I alpha Rabbit pAb at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunocytochemistry/ Immunofluorescence: Protein Kinase A regulatory subunit I alpha Antibody - Azide and BSA Free [NBP2-93193] - Immunofluorescence analysis of NIH/3T3 cells using Protein Kinase A regulatory subunit I alpha Rabbit pAb at dilution of 1:100. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western Blot: Protein Kinase A regulatory subunit I alpha Antibody - Azide and BSA Free [NBP2-93193] - Western blot analysis of lysates from wild type (WT) and Protein Kinase A regulatory subunit I alpha knockout (KO) 293T cells, using [KO Validated] Protein Kinase A regulatory subunit I alpha Rabbit pAb at 1:500 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-93193-0.1ml

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-93193

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

