

Product Datasheet

SARS-CoV-2 Nucleocapsid Antibody (B3451M) - BSA Free NBP3-00513-200ug

Unit Size: 200 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-00513

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-00513



NBP3-00513-200ug

SARS-CoV-2 Nucleocapsid Antibody (B3451M) - BSA Free

Product Information	
Unit Size	200 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	B3451M
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	Phosphate Buffered Saline, pH 7.4
Product Description	
Description	Novus Biologicals Mouse SARS-CoV-2 Nucleocapsid Antibody (B3451M) - BSA Free (NBP3-00513) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	43740575
Gene Symbol	N
Species	SARS-CoV-2, SARS-CoV
Reactivity Notes	This antibody, paired with SARS-CoV-2 Nucleocapsid Antibody (B3449M) [NBP3-00514], does not cross-react with other coronaviruses except SARS-CoV.
Specificity/Sensitivity	Specific for SARS-CoV-2 Nucleocapsid Protein.
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein.
Product Application Details	
Applications	Western Blot, ELISA, ELISA Capture (Matched Antibody Pair), Lateral Flow Assay
Recommended Dilutions	Western Blot, ELISA, Lateral Flow Assay, ELISA Capture (Matched Antibody Pair)
Application Notes	This sandwich ELISA capture antibody paired with the sandwich ELISA detection antibody [NBP3-00514] is used for highly specific and sensitive detection of SARS-CoV-2 Nucleocapsid protein.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-00513-200ug

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-00513

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

