

# Product Datasheet

## Fibrillarin Antibody (4A4) - BSA Free NBP3-05533-100ul

Unit Size: 100 ul

Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-05533](http://www.novusbio.com/NBP3-05533)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-05533](http://www.novusbio.com/reviews/destination/NBP3-05533)



**NBP3-05533-100ul**

Fibrillarin Antibody (4A4) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4A4
Preservative	5mM Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	50% PBS, 50% glycerol
Target Molecular Weight	34.5 kDa


Product Description	
Description	Novus Biologicals Mouse Fibrillarin Antibody (4A4) - BSA Free (NBP3-05533) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2091
Gene Symbol	FBL
Species	Human, Mouse, Rat, Canine, Equine, Monkey
Immunogen	Recombinant full length human Fibrillarin sequence expressed in and purified from E. coli.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, Immunohistochemistry 1:1000-1:5000, Immunocytochemistry/ Immunofluorescence 1:1000-1:5000




## Images

**Western Blot: Fibrillarin Antibody [NBP3-05533]** - Western blot analysis of different tissue and cell line lysates using Fibrillarin Antibody at a dilution of 1:2,000, in green: [1] protein standard (red), [2] rat whole brain, [3] mouse whole brain, [4] NIH-3T3 cells, [5] HEK293, [6] HeLa, and [7] SH-SY5Y cells. Strong band at ~35kDa corresponds to Fibrillarin protein seen in all cell line lysates. A much weaker band is seen in tissue lysates since Fibrillarin is more heavily expressed in rapidly dividing cells. The blot simultaneously was probed with rabbit pAb to UCHL1 at a dilution of 1:3,000, in red, revealing the ~25kDa UCHL1 protein in lysates of tissues containing neurons or cells with neuronal properties.



**Immunocytochemistry/Immunofluorescence: Fibrillarin Antibody [NBP3-05533]** - Immunofluorescent analysis of (A) C6 rat glioma cells and (B) HEK293 human embryonic kidney cells stained with Fibrillarin Antibody at a dilution of 1:1,000 in red, in both cases costained with chicken pAb to vimentin at a dilution of 1:10,000 in green. The blue is DAPI staining of nuclear DNA. The Fibrillarin antibody detects Fibrillarin protein localized in nucleoli while the vimentin antibody produces strong staining of cytoplasmic intermediate filaments.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-05533-100ul**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-05533](http://www.novusbio.com/reviews/submit/NBP3-05533)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

