

Product Datasheet

Human, Mouse, Rat VPS35 ELISA Kit (Colorimetric) NBP3-11045

Unit Size: 96 wells

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-11045

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-11045



NBP3-11045


Human, Mouse, Rat VPS35 ELISA Kit (Colorimetric)

Product Information	
Unit Size	96 wells
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C.
Product Description	
Gene ID	55737
Gene Symbol	VPS35
Species	Human, Mouse, Rat
Kit Components	Anti-VPS35 Immunoassay Plate - 1 plate, 5X VPS35 Extraction Reagent - 10mL, Recombinant VPS35 Standard - 2 vials, Standard and Sample Diluent (Red) - 50mL, 10X Wash Buffer Concentrate - 100mL, Anti-VPS35 Biotinylated Antibody Concentrate - 150uL, Anti-VPS35 Biotinylated Antibody Diluent (Green) - 13mL, 100X Poly HRP Concentrate - 150uL, Streptavidin: HRP Diluent (Purple) - 13mL, TMB Substrate - 13mL, Stop Solution - 13mL
Standard Curve Range	78.125 - 5000 pg/mL
Sensitivity	50 pg/mL
Assay Type	Competitive ELISA
Spiking Recovery	95.5%
Suitable Sample Type	Cell Lysate, Tissue
Sample Volume	100 uL
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA




Images


ELISA: Human, Mouse, Rat VPS35 ELISA Kit (Colorimetric) [NBP3-11045] - Standard Reference Curve



ELISA: Human, Mouse, Rat VPS35 ELISA Kit (Colorimetric) [NBP3-11045] - Linearity was determined by using low (rat kidney) and high (mouse kidney) natural samples, diluted 1/100 and mixed in various proportions.



ELISA: Human, Mouse, Rat VPS35 ELISA Kit (Colorimetric) [NBP3-11045] - VPS35 Determination of WT-A549 vs KN-549 Lysates Graph for the VPS35 Kit. Assay Type: Competitive ELISA. Detection Method: Colorimetric Assay.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-11045

H00055737-P01-10ug	Recombinant Human VPS35 GST (N-Term) Protein
NBP2-13075	beta Amyloid Antibody (MOAB-2) - BSA Free
NB100-1397	VPS35 Antibody
NB300-268	LRRK2 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-11045

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

