

# Product Datasheet

## FCRL3/FcRH3 Antibody (H5) - BSA Free NBP3-11831

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-11831](http://www.novusbio.com/NBP3-11831)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-11831](http://www.novusbio.com/reviews/destination/NBP3-11831)



**NBP3-11831**

FCRL3/FcRH3 Antibody (H5) - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	H5
Preservative	15mM Sodium Azide
Isotype	IgG2b Kappa
Purity	Protein A purified
Buffer	Phosphate buffered saline (PBS), pH 7.4

Product Description	
Description	Novus Biologicals Mouse FCRL3/FcRH3 Antibody (H5) - BSA Free (NBP3-11831) is a monoclonal antibody validated for use in Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	115352
Gene Symbol	FCRL3
Species	Human
Specificity/Sensitivity	The mouse monoclonal antibody H5 recognizes an epitope within extracellular part of FCRL3/FcRH3, a transmembrane glycoprotein expressed mainly on NK cells, and T and B cell subsets.
Immunogen	DNA-immunization followed by a boost with FCRL3/FcRH3 transfected cells

Product Application Details	
Applications	Flow Cytometry
Recommended Dilutions	Flow Cytometry 1-4 ug/ml

**Images**

Flow Cytometry: FCRL3/FcRH3 Antibody (H5) [NBP3-11831] - Multicolor surface staining of human lymphocytes stained using anti-human CD19 (LT19) APC antibody (10 ul reagent / 100 ul of peripheral whole blood) and anti-human CD307c (H5) purified antibody (1,7 ug/ml, GAM-FITC).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-11831**

---

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP2-62615PEP	FCRL3/FcRH3 Recombinant Protein Antigen

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-11831](http://www.novusbio.com/reviews/submit/NBP3-11831)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

