

# Product Datasheet

## Cytokeratin 7 Antibody (OV-TL 12/30) [DyLight 350] NBP3-12013UV

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-12013UV](http://www.novusbio.com/NBP3-12013UV)

Updated 5/1/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-12013UV](http://www.novusbio.com/reviews/destination/NBP3-12013UV)



**NBP3-12013UV**

Cytokeratin 7 Antibody (OV-TL 12/30) [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                    |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                       |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                     |
| Clone                       | OV-TL 12/30                                                                                                                                                                                                                                                                                                    |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                             |
| Isotype                     | IgG1 Kappa                                                                                                                                                                                                                                                                                                     |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                    |
| Purity                      | Protein A purified                                                                                                                                                                                                                                                                                             |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                             |
| Product Description         |                                                                                                                                                                                                                                                                                                                |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                          |
| Gene ID                     | 3855                                                                                                                                                                                                                                                                                                           |
| Gene Symbol                 | KRT7                                                                                                                                                                                                                                                                                                           |
| Species                     | Human                                                                                                                                                                                                                                                                                                          |
| Specificity/Sensitivity     | This antibody recognizes cytokeratin-7. Cytokeratins are a subfamily of intermediate filament proteins but this antibody does not recognize other intermediate filament proteins. The antibody reacts with proteins that are found in most ductal, glandular, transitional, and biliary duct epithelial cells. |
| Immunogen                   | This antibody was originally raised by immunizing BALB/c mice with ovarian carcinoma cell line OTN11.                                                                                                                                                                                                          |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                              |
| Product Application Details |                                                                                                                                                                                                                                                                                                                |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                                                                                                                                                                                                    |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence                                                                                                                                                                                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-12013UV**

---

|                    |                                                             |
|--------------------|-------------------------------------------------------------|
| NBP1-43319UV-0.5ml | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 350] |
| H00003855-P01-10ug | Recombinant Human Cytokeratin 7 GST (N-Term) Protein        |
| 8184-CK-050        | Choline Kinase beta [Unconjugated]                          |
| NBL1-12396         | Cytokeratin 7 Overexpression Lysate                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-12013UV](http://www.novusbio.com/reviews/submit/NBP3-12013UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

