

# Product Datasheet

## Common gamma Chain/IL-2 R gamma Antibody (3E12) [PE] NBP3-12039PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-12039PE](http://www.novusbio.com/NBP3-12039PE)

Updated 5/1/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-12039PE](http://www.novusbio.com/reviews/destination/NBP3-12039PE)



**NBP3-12039PE**

Common gamma Chain/IL-2 R gamma Antibody (3E12) [PE]

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                              |
|------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                                                                       |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                  |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                                                                                     |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                                                                                                   |
| <b>Clone</b>                       | 3E12                                                                                                                                                                                                                                                                         |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                                                                           |
| <b>Isotype</b>                     | IgG2b Kappa                                                                                                                                                                                                                                                                  |
| <b>Conjugate</b>                   | PE                                                                                                                                                                                                                                                                           |
| <b>Purity</b>                      | Protein A purified                                                                                                                                                                                                                                                           |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                                                                                                          |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                              |
| <b>Host</b>                        | Rat                                                                                                                                                                                                                                                                          |
| <b>Gene ID</b>                     | 3561                                                                                                                                                                                                                                                                         |
| <b>Gene Symbol</b>                 | IL2RG                                                                                                                                                                                                                                                                        |
| <b>Species</b>                     | Mouse                                                                                                                                                                                                                                                                        |
| <b>Specificity/Sensitivity</b>     | This conformation-dependent antibody binds to the murine common gamma-chain, shared by the IL-2, IL-4, IL-7, IL-9, IL-15, and IL-21 receptor complexes. It recognises an epitope located in a membrane-proximal portion of the common gamma-chain, near Arg-197 of domain 2. |
| <b>Immunogen</b>                   | This antibody was raised by immunising Lewis rats with transfected rat myeloma cell line YB2/0, expressing murine cytoplasmic-tailless IL-2R common gamma-chain (YB2/0-myt).                                                                                                 |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                              |
| <b>Applications</b>                | Flow Cytometry, Immunoprecipitation, Block/Neutralize                                                                                                                                                                                                                        |
| <b>Recommended Dilutions</b>       | Flow Cytometry, Immunoprecipitation, Block/Neutralize                                                                                                                                                                                                                        |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                       |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-12039PE**

---

|                |                                                             |
|----------------|-------------------------------------------------------------|
| NBP1-43323PE   | Rat IgG2b Kappa Light Chain Isotype Control (149/10H5) [PE] |
| 6507-IL-010/CF | IL-4 [Unconjugated]                                         |
| 384-RG-050     | Common gamma Chain/IL-2 R gamma [Unconjugated]              |
| 202-IL-010     | IL-2 [Unconjugated]                                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-12039PE](http://www.novusbio.com/reviews/submit/NBP3-12039PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

