

Product Datasheet

SARS-CoV-2 Nucleocapsid Antibody (019) [DyLight 350] NBP3-12735UV

Unit Size: 100 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12735UV

Updated 7/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12735UV



NBP3-12735UV**SARS-CoV-2 Nucleocapsid Antibody (019) [DyLight 350]**

| Product Information | |
|------------------------------------|--|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 019 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 350 |
| Purity | Protein A purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rabbit |
| Gene ID | 43740575 |
| Gene Symbol | N |
| Species | SARS-CoV |
| Reactivity Notes | No cross-reactivity in ELISA with: MERS-CoV Nucleoprotein protein HCoV-229E Nucleoprotein protein HCoV-NL63 Nucleoprotein protein HCoV-HKU1(isolate N5) Nucleoprotein protein HCoV-OC43 Nucleoprotein |
| Specificity/Sensitivity | Has cross-reactivity in ELISA with: SARS-CoV-2 Omicron (B.1.1.529) Nucleocapsid Protein No cross-reactivity in ELISA with: MERS-CoV Nucleocapsid protein HCoV-229E Nucleocapsid protein HCoV-NL63 Nucleocapsid protein HCoV-HKU1(isolate N5) Nucleocapsid protein HCoV-OC43 Nucleocapsid protein |
| Immunogen | Recombinant SARS-CoV-2 Nucleocapsid Protein |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, ELISA |
| Recommended Dilutions | Western Blot, ELISA |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-12735UV

| | |
|--------------|---|
| NBP2-24891UV | Rabbit IgG Isotype Control [DyLight 350] |
| NBP2-90975 | Recombinant SARS-CoV-2 Nucleocapsid His (C-Term) Protein |
| NBP3-49916 | SARS-CoV-2 SARS-CoV-2 Nucleocapsid ELISA Kit (Colorimetric) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12735UV

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

