

Product Datasheet

Beta Dystroglycan Antibody (GT9112) - Azide and BSA Free NBP3-13642

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13642

Updated 3/3/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13642



NBP3-13642

Beta Dystroglycan Antibody (GT9112) - Azide and BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GT9112
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS


Product Description	
Description	Centrifuge briefly prior to opening.
Host	Mouse
Gene ID	1605
Gene Symbol	DAG1
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human Beta Dystroglycan / beta Dystroglycan. The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin




Images


Western Blot: Beta Dystroglycan Antibody (GT9112) [NBP3-13642] - Wild-type (WT) and Beta Dystroglycan knockout (KO) 293T cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with Beta Dystroglycan antibody [GT9112] (NBP3-13642) diluted at 1:500. The HRP-conjugated anti-mouse IgG antibody (NBP2-19382) was used to detect the primary antibody.



Immunohistochemistry-Paraffin: Beta Dystroglycan Antibody (GT9112) [NBP3-13642] - Beta Dystroglycan antibody [GT9112] detects Beta Dystroglycan protein at cell membrane and cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse placenta. Beta Dystroglycan stained by Beta Dystroglycan antibody [GT9112] (NBP3-13642) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: Beta Dystroglycan Antibody (GT9112) [NBP3-13642] - Various whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with Beta Dystroglycan antibody [GT9112] (NBP3-13642) diluted at 1:500. The HRP-conjugated anti-mouse IgG antibody (NBP2-19382) was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13642

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
236-EG-200	EGF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13642

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

