

Product Datasheet

CA19-9/Sialyl Lewis A Antibody (GT252) - Azide and BSA Free NBP3-13651

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13651

Updated 9/25/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13651



NBP3-13651

CA19-9/Sialyl Lewis A Antibody (GT252) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	GT252
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS, 20% Glycerol

Product Description	
Description	Centrifuge briefly prior to opening.
Host	Mouse
Species	Human
Immunogen	The immunogen used to generate this antibody corresponds to CA19-9/Sialyl Lewis A.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Sandwich ELISA
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Sandwich ELISA

Images


Western Blot: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membranes were blotted with CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) diluted at 1:1000 and competitor's antibody diluted at 1:250. The HRP-conjugated anti-mouse IgG antibody (NBP2-19382) was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - CA19-9/Sialyl Lewis A antibody [GT252] detects CA19-9/Sialyl Lewis A protein at cell membrane by immunofluorescent analysis. Sample: SW1116 (high expression) and HT29 (low expression) cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: CA19-9/Sialyl Lewis A stained by CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) diluted at 1:1000 and competitor's antibody diluted at 1:1000. Blue: Fluoroshield with DAPI.




Flow Cytometry: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) detects CA19-9/Sialyl Lewis A by flow cytometry analysis. Sample: COLO 205 cell. Black: Unlabelled sample was used as a control. Red: CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) dilution: 1:50. Blue: Highly cited antibody [116-NS-19-9] dilution: 1:50. Acquisition of 20,000 events were collected for FACS analysis.




Western Blot: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) diluted at 1:1000. The HRP-conjugated anti-mouse IgG antibody (NBP2-19382) was used to detect the primary antibody.




ELISA: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - Competition analysis of CA 19-9 antibody [GT252] (NBP3-13651) (330 ng/mL) binding to immobilized recombinant HSA-CA 19-9 (coated at 3 ug/mL) in the presence of soluble Lewis A or Lewis X tetrasaccharides (1.0-0.01 mM). Bound antibody was detected by Goat Anti-Mouse IgG antibody (HRP) (NBP2-19382) (1:10000).



Sandwich ELISA: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - Sandwich ELISA detection of Colo 205 and Jurkat whole cell extracts using CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) as capture antibody at concentration of 5 ug/mL and HRP-conjugated Ca19-9 antibody [GT8510] (NBP3-13652) as detection antibody at concentration of 1 ug/mL. Please notice that NBP3-13652 needs to be conjugated to HRP to function as the detection antibody when paired with NBP3-13651. Please contact us for custom HRP-conjugated antibody.



Immunohistochemistry-Paraffin: CA19-9/Sialyl Lewis A Antibody (GT252) [NBP3-13651] - CA19-9/Sialyl Lewis A antibody [GT252] detects CA19-9/Sialyl Lewis A protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded human pancreatic cancer. CA19-9/Sialyl Lewis A stained by CA19-9/Sialyl Lewis A antibody [GT252] (NBP3-13651) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13651

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-59963-1Kit	Human CA19-9/Sialyl Lewis A ELISA Kit (Colorimetric)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13651

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

