

# Product Datasheet

## FoxP3 Antibody (BLR034F) NBP3-14670

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14670](http://www.novusbio.com/NBP3-14670)

Updated 2/24/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14670](http://www.novusbio.com/reviews/destination/NBP3-14670)



**NBP3-14670**

FoxP3 Antibody (BLR034F)


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 2 - 8 C / 1 year from date of receipt
Clonality	Monoclonal
Clone	BLR034F
Preservative	0.09% Sodium Azide
Purity	>95%
Buffer	Borate Buffered Saline (BBS) pH 8.2 with , BSA Free

Product Description	
Description	Novus Biologicals Rabbit FoxP3 Antibody (BLR034F) (NBP3-14670) is a recombinant monoclonal antibody validated for use in Multiplex Immunofluorescence, IHC, WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	50943
Gene Symbol	FOXP3
Species	Human
Immunogen	between 275 and 325


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Multiplex Immunofluorescence
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100 - 1:500, Immunoprecipitation 20 ul/1mg lysate, Immunohistochemistry-Paraffin 1:100 - 1:500, Multiplex Immunofluorescence 1:250
Application Notes	Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections.

**Images**


Western Blot: FoxP3 Antibody (BLR034F) [NBP3-14670] - Whole cell lysate (50 ug) from Jurkat, HEK293T, MJ, MOLT-4, and RPMI-8226 cells prepared using NETN lysis buffer. Rabbit anti-FOXP3 recombinant monoclonal antibody (NBP3-14670 lot 2) used at 1:1000. Secondary: HRP-conjugated goat anti-rabbit IgG (NB7156). Chemiluminescence with an exposure time of 30 seconds. Lower Panel: Rabbit anti-COPB2 antibody.




Immunohistochemistry: FoxP3 Antibody (BLR034F) [NBP3-14670] - FFPE section of human breast carcinoma. Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody used at 1:250. Secondary: HRP-conjugated goat antirabbit IgG. Substrate: Opal™. Counterstain: DAPI (blue).




Immunohistochemistry-Paraffin: FoxP3 Antibody (BLR034F) [NBP3-14670] - Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.




Immunohistochemistry-Paraffin: FoxP3 Antibody (BLR034F) [NBP3-14670] - Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.



Immunohistochemistry: FoxP3 Antibody (BLR034F) [NBP3-14670] - Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.



Immunoprecipitation: FoxP3 Antibody (BLR034F) [NBP3-14670] -  
Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded)  
from MJ cells prepared using NETN lysis buffer. Antibodies: Rabbit anti-  
FOXP3 recombinant monoclonal antibody used for IP at 20  $\mu$ l/mg lysate.  
FOXP3 was also immunoprecipitated by a previous lot of this antibody.  
For blotting immunoprecipitated FOXP3, was used at 1:1000.  
Chemiluminescence with an exposure time of 30 seconds.





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to NBP3-14670

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NB600-246PEP	FoxP3 Antibody Blocking Peptide

---

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14670](http://www.novusbio.com/reviews/submit/NBP3-14670)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

