

Product Datasheet

eIF4E Antibody (7A7P2)

NBP3-15629-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15629

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15629



NBP3-15629-100ul

eIF4E Antibody (7A7P2)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	7A7P2
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Target Molecular Weight	25 kDa


Product Description	
Description	Novus Biologicals Rabbit eIF4E Antibody (7A7P2) (NBP3-15629) is a recombinant monoclonal antibody validated for use in WB, ELISA and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1977
Gene Symbol	EIF4E
Species	Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 118-217 of human eIF4E (P06730). NKQQRSDLDLDRFWLETLLCLIGESFDDYSDDVCGAVVNVRAKGDKIAIWTTEC ENREAVTHIGRVYKERLGLPPKIVIGYQSHADTATKSGSTTKNRFVV

Product Application Details	
Applications	Western Blot, ELISA, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500 - 1:1000, ELISA, Immunoprecipitation 0.5µg-4µg antibody for 400µg-600µg extracts of whole cells




Images

Western Blot: eIF4E Antibody (7A7P2) [NBP3-15629] - Western blot analysis of extracts of various cell lines, using eIF4E antibody (NBP3-15629) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 3min.



Immunoprecipitation: eIF4E Antibody (7A7P2) [NBP3-15629] - Immunoprecipitation analysis of 600ug extracts of Mouse testis using 3ug eIF4E antibody (NBP3-15629). Western blot was performed from the immunoprecipitate using eIF4E (NBP3-15629) at a dilution of 1:500.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15629-100ul

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15629

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

