

Product Datasheet

Complement C3a Antibody (6E8) - BSA Free NBP3-17998-100ul

Unit Size: 100 ul

Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-17998

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-17998



NBP3-17998-100ul

Complement C3a Antibody (6E8) - BSA Free


Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6E8
Preservative	0.035% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	50% PBS, 50% glycerol
Target Molecular Weight	185 kDa

Product Description	
Description	Novus Biologicals Mouse Complement C3a Antibody (6E8) - BSA Free (NBP3-17998) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	718
Gene Symbol	C3
Species	Human
Immunogen	Recombinant netrin domain of human C3. UniProt: P01024

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:1000 - 1:5000, ELISA

Images

Western Blot: Complement C3a Antibody (6E8) [NBP3-17998] - Western blot analysis of Complement C3a antibody (6E8) on human serum samples. 0.1ug purified human C3 protein (lane 1) and 10 ug normal human serum protein (lane 2) were blotted with Complement C3a antibody (6E8) at a concentration 1:3000. The Complement C3a antibody (6E8) monoclonal binds strongly and cleanly to a band at about 115 kDa which represents the intact alpha subunit of C3 and a weaker proteolytic band at approximately 40kDa. Bands at 190kDa and above are likely the pro-C3 and its glycosylated form.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-17998-100ul

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-17998

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

