

# Product Datasheet

## metallothionein Antibody (2B5) - BSA Free NBP3-18262

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-18262](http://www.novusbio.com/NBP3-18262)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-18262](http://www.novusbio.com/reviews/destination/NBP3-18262)



**NBP3-18262**

metallothionein Antibody (2B5) - BSA Free


Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2B5
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS pH 7.4, 50% glycerol

Product Description	
Description	Novus Biologicals Mouse metallothionein Antibody (2B5) - BSA Free (NBP3-18262) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Species	Bacteria
Reactivity Notes	Bacteria ( <i>Pseudomonas aeruginosa</i> )
Specificity/Sensitivity	Detects 36 kDa band corresponding to <i>Pseudomonas aeruginosa</i> Metallothionein.
Immunogen	Recombinant full length <i>Pseudomonas aeruginosa</i> Metallothionein protein

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:1000, ELISA 1:800
Application Notes	Optimal dilutions for assays should be determined by the user.

**Images**

Western Blot: metallothionein Antibody (2B5) [NBP3-18262] - Western Blot analysis of *Pseudomonas aeruginosa* metallothionein (PmtA) GST tagged showing detection of 36 kDa metallothionein protein using Mouse Anti-metallothionein Monoclonal Antibody, Clone 2B5 (NBP3-18262). Lane 1: Molecular Weight Ladder (MW). Lane 2: GST. Lane 3: Purified *Pseudomonas aeruginosa* metallothionein (PmtA) GST tagged. Load: 1 ug. Block: 5% Skim Milk powder in TBST. Primary Antibody: Mouse Anti-metallothionein Monoclonal Antibody (NBP3-18262) at 1:1000 for 2 hours at RT. Secondary Antibody: Goat Anti-Mouse IgG:HRP at 1:5000 for 1 hour at RT. Color Development: ECL solution for 5 min in RT. Predicted/Observed Size: 36 kDa.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-18262**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-18262](http://www.novusbio.com/reviews/submit/NBP3-18262)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

