

Product Datasheet

Integrin alpha 5/CD49e Antibody (SAM1) - BSA Free NBP3-18529

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18529

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18529




NBP3-18529**Integrin alpha 5/CD49e Antibody (SAM1) - BSA Free**


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	SAM1
Preservative	15mM Sodium Azide
Isotype	IgG2b
Purity	Protein A purified
Buffer	Phosphate buffered saline (PBS), pH 7.4
Product Description	
Description	Novus Biologicals Mouse Integrin alpha 5/CD49e Antibody (SAM1) - BSA Free (NBP3-18529) is a monoclonal antibody validated for use in IHC, Flow and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3678
Gene Symbol	ITGA5
Species	Human, Primate
Immunogen	U937 cells
Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Flow Cytometry 1-4 ug/ml, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Frozen

Images

Flow Cytometry: Integrin alpha 5/CD49e Antibody (SAM1) [NBP3-18529]
- Separation of human Integrin alpha 5/CD49e positive monocytes (red-filled) from human Integrin alpha 5/CD49e negative lymphocytes (black-dashed) in flow cytometry analysis (surface staining) of peripheral whole blood stained using anti-human Integrin alpha 5/CD49e (SAM1) purified antibody (concentration in sample 1.7 ug/ml, GAM APC).



Flow Cytometry: Integrin alpha 5/CD49e Antibody (SAM1) [NBP3-18529]
- Flow cytometry surface staining pattern of human peripheral whole blood stained using anti-human Integrin alpha 5/CD49e (SAM1) purified antibody (concentration in sample 1.7 ug/ml, GAM APC).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-18529

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)
NBP1-84576PEP	Integrin alpha 5/CD49e Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18529

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

