

Product Datasheet

LSD1 Antibody - BSA Free NBP3-18684

Unit Size: 50 ug

Store at -20C short term. Aliquot and store at -80C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18684

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18684



NBP3-18684

LSD1 Antibody - BSA Free


Product Information	
Unit Size	50 ug
Concentration	2.2 mg/ml
Storage	Store at -20C short term. Aliquot and store at -80C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide and 0.05% ProClin 300
Purity	Affinity purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit LSD1 Antibody - BSA Free (NBP3-18684) is a polyclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and CHIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	23028
Gene Symbol	KDM1A
Species	Human
Immunogen	Polyclonal antibody raised in rabbit against human LSD1 (Lysine-specific demethylase 1), using a KLH-conjugated synthetic peptide from the inner part of the protein.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Chromatin Immunoprecipitation (ChIP), Chromatin Immunoprecipitation Sequencing
Recommended Dilutions	Western Blot 1:4000, ELISA 1:1000-1:10000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:200, Immunohistochemistry-Paraffin 1:200-1:500, Chromatin Immunoprecipitation (ChIP) 1 ug/ChIP, Chromatin Immunoprecipitation Sequencing 1 ug/ChIP
Application Notes	Please note that the optimal antibody amount per IP should be determined by the end-user. We recommend testing 1-5 ug per IP.

Images


Western Blot: LSD1 Antibody [NBP3-18684] - Figure 5. Western blot analysis using the antibody directed against LSD1. Whole cell extracts (40 ug) from HeLa cells transfected with LSD1 siRNA (lane 2) and from an untransfected control (lane 1) were analysed by Western blot using the antibody against LSD1 diluted 1:5,000 in TBS-Tween containing 5% skimmed milk. The position of the protein of interest is indicated on the right; the marker (in kDa) is shown on the left.




Immunocytochemistry/Immunofluorescence: LSD1 Antibody [NBP3-18684] - Figure 6. Immunofluorescence using the antibody directed against LSD1 HeLa cells were stained with the antibody against LSD1 and with DAPI. Cells were fixed with 4% formaldehyde for 10 and blocked with PBS/TX-100 containing 5% normal goat serum and 1% BSA. The cells were immunofluorescently labelled with the LSD1 antibody (left) diluted 1:200 in blocking solution followed by an anti-rabbit antibody conjugated to Alexa488. The middle panel shows staining of the nuclei with DAPI. A merge of the two stainings is shown on the right.




Chromatin Immunoprecipitation Sequencing: LSD1 Antibody [NBP3-18684] - Figure 2. ChIP-seq results obtained with the antibody directed against LSD1 ChIP was performed on sheared chromatin from 4,000,000 K562 cells using 1 ug of the antibody against LSD1 as described above. The IPd DNA was subsequently analysed on an Illumina HiSeq. Library preparation, cluster generation and sequencing were performed according to the manufacturers instructions. The 50 bp tags were aligned to the human genome using the BWA algorithm. Figure 2 shows the peak distribution along the complete sequence and a 600 kb region of the X-chromosome (figure 2A and B) and in three regions surrounding the MYT1, RBM19 and TGFBR3 positive control genes, respectively (figure 2C, D and E). The position of the amplicon used for ChIP-qPCR is indicated by an arrow.




Western Blot: LSD1 Antibody [NBP3-18684] - Figure 4. Western blot analysis using the antibody directed against LSD1 Western blot was performed using nuclear extracts from HeLa cells (40 ug) and the antibody against LSD1 diluted 1:4,000 in TBS- Tween containing 5% skimmed milk. The molecular weight marker (in kDa) is shown on the left. The location of the protein of interest is indicated on the right.



Chromatin Immunoprecipitation Sequencing: LSD1 Antibody [NBP3-18684] - Figure 1. ChIP results obtained with the antibody directed against LSD1 ChIP was performed with the antibody against LSD1 on sheared chromatin from 4,000,000 K562 cells using the "iDeal ChIP-seq" kit .. An antibody titration consisting of 1, 2, 5 and 10 ug per ChIP experiment was analysed. IgG (2 ug/IP) was used as negative IP control. QPCR was performed with primers for specific regions in the MYT1, RBM19, and TGFBR3 genes, used as positive controls, and for the MYOD1 gene, used as negative control. Figure 1 shows the recovery, expressed as a % of input (the relative amount of immunoprecipitated DNA compared to input DNA after qPCR analysis).



ELISA: LSD1 Antibody [NBP3-18684] - Figure 3. Determination of the antibody titer To determine the titer of the antibody, an ELISA was performed using a serial dilution of the antibody directed against LSD1 in antigen coated wells. By plotting the absorbance against the antibody dilution (Figure 3), the titer of the antibody was estimated to be 1:176,000.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-18684

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
H00023028-P01-10ug	Recombinant Human LSD1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18684

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

