

Product Datasheet

UBA3 Antibody (S05-6C8)

NBP3-20008-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20008

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20008



NBP3-20008-100ul

UBA3 Antibody (S05-6C8)


| Product Information | |
|-------------------------|--|
| Unit Size | 100 ul |
| Concentration | 0.3 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | S05-6C8 |
| Preservative | 0.01% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA |
| Target Molecular Weight | 52 kDa |

| Product Description | |
|---------------------|---|
| Description | Novus Biologicals Rabbit UBA3 Antibody (S05-6C8) (NBP3-20008) is a recombinant monoclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 9039 |
| Gene Symbol | UBA3 |
| Species | Human, Mouse, Rat |
| Immunogen | A synthetic peptide of human UBA3 (Uniprot # Q8TBC4) |

| Product Application Details | |
|-----------------------------|-----------------------------|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot 1/2000-1/10000 |

Images

Western Blot: UBA3 Antibody (S05-6C8) [NBP3-20008] - Western blot detection of UBA3 in HeLa, A549, HL-60 cell lysates using NBP3-20008 (1:1000 diluted). Predicted band size: 52kDa. Observed band size: 52kDa.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-20008-100ul

| | |
|-------------|---|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20008

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

