

Product Datasheet

KIR3DL3/CD158z Antibody (CH21) - BSA Free NBP3-21499

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-21499

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-21499



NBP3-21499**KIR3DL3/CD158z Antibody (CH21) - BSA Free**

| Product Information | |
|----------------------------|---|
| Unit Size | 0.1 mg |
| Concentration | 1 mg/ml |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Monoclonal |
| Clone | CH21 |
| Preservative | 15mM Sodium Azide |
| Isotype | IgG2a |
| Purity | Protein A purified |
| Buffer | Phosphate buffered saline (PBS), pH 7.4 |


| Product Description | |
|----------------------------|---|
| Description | Novus Biologicals Mouse KIR3DL3/CD158z Antibody (CH21) - BSA Free (NBP3-21499) is a monoclonal antibody validated for use in Flow and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Mouse |
| Gene ID | 115653 |
| Gene Symbol | KIR3DL3 |
| Species | Human |
| Immunogen | Human KIR3DL3/CD158z transfectants |

| Product Application Details | |
|------------------------------------|--|
| Applications | Flow Cytometry, Immunoprecipitation |
| Recommended Dilutions | Flow Cytometry 3-12 ug/ml, Immunoprecipitation |




Images

Flow cytometry surface staining pattern of KIR3DL3/CD158z transfected HEK-293 cells co-transfected with YFP coding plasmid using anti-human KIR3DL3/CD158z (CH21) purified antibody (concentration in sample 10 ug/ml, GAM APC).



KIR3DL3/CD158z purified / GAM APC

Flow cytometry surface staining patterns of non-transfected HEK-293 cells and HEK-293 cells transfected with KIR-family coding plasmids co-transfected with YFP coding plasmid using anti-human KIR3DL3/CD158z (CH21) purified antibody (concentration in sample 10 ug/ml, GAM APC).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-21499

| | |
|-------------|--|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB7539 | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-96778 | Mouse IgG2a Isotype Control (M2A) |
| 3015-KR-050 | KIR2DL2/CD158b1 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-21499

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

