

# Product Datasheet

## uPAR Antibody (HL2169) - Azide and BSA Free NBP3-25734

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-25734](http://www.novusbio.com/NBP3-25734)

Updated 3/4/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-25734](http://www.novusbio.com/reviews/destination/NBP3-25734)



**NBP3-25734**

uPAR Antibody (HL2169) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2169
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS
Target Molecular Weight	45 kDa

Product Description	
Host	Rabbit
Gene ID	5329
Gene Symbol	PLAUR
Species	Human, Mouse, Rat
Immunogen	Recombinant fragment of human uPAR.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry Assay dependent, Immunocytochemistry/ Immunofluorescence Assay dependent, Immunohistochemistry-Paraffin Assay dependent

**Images**


Immunohistochemistry-Paraffin: uPAR Antibody (HL2169) - Azide and BSA Free [NBP3-25734] - uPAR antibody [HL2169] detects uPAR protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissues. uPAR stained by uPAR antibody [HL2169] (NBP3-25734) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Western Blot: uPAR Antibody (HL2169) - Azide and BSA Free [NBP3-25734] - Wild-type (WT) and PLAUR knockout (KO) HeLa cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with uPAR antibody [HL2169] (NBP3-25734) diluted at 1:250. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.




Western Blot: uPAR Antibody (HL2169) - Azide and BSA Free [NBP3-25734] - Untreated (-) and treated (+) U937 whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with uPAR antibody [HL2169] (NBP3-25734) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.



Immunohistochemistry-Paraffin: uPAR Antibody (HL2169) - Azide and BSA Free [NBP3-25734] - uPAR antibody [HL2169] detects uPAR protein by immunohistochemical analysis. Sample: Paraffin-embedded rat tissues. uPAR stained by uPAR antibody [HL2169] (NBP3-25734) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunocytochemistry/Immunofluorescence: uPAR Antibody (HL2169) - Azide and BSA Free [NBP3-25734] - uPAR antibody [HL2169] detects uPAR protein at cell membrane and cytoplasmic vesicles by immunofluorescent analysis. Sample: U937 cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: uPAR stained by uPAR antibody [HL2169] (NBP3-25734) diluted at 1:500. Blue: Fluoroshield with DAPI .





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-25734**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-52381-0.05mg	Recombinant Human uPAR His Protein

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-25734](http://www.novusbio.com/reviews/submit/NBP3-25734)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

