

Product Datasheet

Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free NBP3-25995

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25995

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25995



NBP3-25995

Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free


Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	RM166
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Rabbit Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free (NBP2-61538) is a recombinant monoclonal antibody validated for use in Multiplex Immunofluorescence, IHC, WB, ELISA, ICC/IF and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	126961
Gene Symbol	H3C14
Species	Human, Vertebrate
Specificity/Sensitivity	This Histone H3 [ac Lys18] antibody (RM166) reacts to Histone H3 acetylated at Lysine 18. No cross reactivity with other acetylated Lysines in Histone H3.
Immunogen	This Histone H3 [ac Lys18] antibody (RM166) was raised against an acetyl-peptide corresponding to the Histone H3 [ac Lys18]


Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Chromatin Immunoprecipitation (ChIP), Multiplex Immunofluorescence
Recommended Dilutions	Western Blot 0.5 ug/ml-2 ug/ml, ELISA 0.2 ug/mL - 1 ug/mL, Immunohistochemistry 1 ug/mL-10 ug/mL, Immunocytochemistry/ Immunofluorescence 0.5 ug/ml - 2 ug/ml, Chromatin Immunoprecipitation (ChIP) 2 ug/mL-10 ug/mL, Multiplex Immunofluorescence

Images


Immunohistochemistry-Paraffin: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - Immunohistochemistry staining of HepG2 cells using Clone RM166.




Western Blot: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - Western Blot of acid extracts from HeLa cells untreated (-) or treated (+) with sodium butyrate, using RM166 at 0.5 ug/mL.




Multiplex Immunoassay: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - RM166 specifically reacts to Histone H3 acetylated at Lysine 18 (K18ac). No cross reactivity with acetylated Lysine 4 (K4ac), Lysine 9 (K9ac), Lysine 14 (K14ac), Lysine 23 (K23ac), Lysine 27 (K27ac), Lysine 36 (K36ac), Lysine 56 (K56ac), lysine 79 (K79ac), or Lysine 122 (K122) in histone H3.



Sandwich ELISA: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - Sandwich ELISA against acetylated histone H3 at Lys 18 using HeLa whole cell lysate, treated or untreated with Sodium Butyrate. Using an anti-Histone H3 (1 ug/mL) as the capture antibody and biotinylated clone RM166 (2 ug/mL) as the detection antibody.




Chromatin Immunoprecipitation: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - ChIP performed on HeLa cells with or without Sodium Butyrate treatment, using Clone RM166 (5ug). Real-time PCR was performed using primers specific to the gene indicated.



Immunocytochemistry/ Immunofluorescence: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - Immunocytochemistry of HeLa cells treated with sodium butyrate, using Clone RM166 (red). Actin filaments have been labeled with fluorescein phalloidin (green).



Chromatin Immunoprecipitation: Histone H3 [ac Lys18] Antibody (RM166) - Azide and BSA Free [NBP3-25995] - ChIP performed on HeLa cells with or without Sodium Butyrate treatment, using Clone RM166 (5ug). Real-time PCR was performed using primers specific to the gene indicated.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25995

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25995

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

