

Product Datasheet

Histone H3 [Monomethyl Lys56] Antibody (RM180) - Azide and BSA Free NBP3-26023

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26023

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26023



NBP3-26023

Histone H3 [Monomethyl Lys56] Antibody (RM180) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	RM180
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit Histone H3 [Monomethyl Lys56] Antibody (RM180) - Azide and BSA Free (NBP2-61529) is a recombinant monoclonal antibody validated for use in Multiplex Immunofluorescence, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	126961
Gene Symbol	H3C14
Species	Human, Vertebrate
Specificity/Sensitivity	This Histone H3 [Monomethyl Lys56] antibody (RM180) reacts to Histone H3 monomethylated at Lysine 56. No cross reactivity with nonmodified Lysine 56 or other methylations in histone H3.
Immunogen	This Histone H3 [Monomethyl Lys56] antibody (RM180) was raised against a monomethyl-peptide corresponding to Histone H3 [Monomethyl Lys56]

Product Application Details	
Applications	Western Blot, ELISA, Multiplex Immunofluorescence
Recommended Dilutions	Western Blot 1 ug/ml - 2 ug/ml, ELISA 0.2 ug/mL - 1 ug/mL, Multiplex Immunofluorescence




Images

Multiplex Immunoassay: Histone H3 [Monomethyl Lys56] Antibody (RM180) - Azide and BSA Free [NBP3-26023] - RM180 specifically reacts to Histone H3 monomethylated at Lysine 56. No cross reactivity with nonmodified Lysine 56, or other methylations in histone H3.



Western Blot: Histone H3 [Monomethyl Lys56] Antibody (RM180) - Azide and BSA Free [NBP3-26023] - Western Blot of acid extracts of HeLa cells (1) and recombinant Histone H3.3 (2), using RM180 at 1 ug/mL, showed a band of Histone H3 monomethylated at Lysine 56 in HeLa cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26023

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26023

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

