

Product Datasheet

5T4 Antibody (13C3) - BSA Free NBP3-26274-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26274

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26274



NBP3-26274-100ul

5T4 Antibody (13C3) - BSA Free

| Product Information | |
|------------------------------------|--|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -20 to -70C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | 13C3 |
| Preservative | 0.03% Proclin 300 |
| Isotype | IgG1 |
| Purity | Affinity purified |
| Buffer | 50% Glycerol, 0.01M PBS, PH 7.4 |
| Product Description | |
| Description | Novus Biologicals Human 5T4 Antibody (13C3) - BSA Free (NBP3-26274) is a recombinant monoclonal antibody validated for use in ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Human |
| Gene ID | 7162 |
| Gene Symbol | TPBG |
| Species | Human, Rhesus Macaque |
| Immunogen | Recombinant Human 5T4 protein [UniProt Q13641] |
| Product Application Details | |
| Applications | ELISA |
| Recommended Dilutions | ELISA |




Images

ELISA: 5T4 Antibody (13C3) [NBP3-26274] - The Binding Activity of 5T4 with 5T4 Antibody.

Activity: Measured by its binding ability in a functional ELISA.


Immobilized Human 5T4 at 2 ug/mL can bind 5T4 Antibody, the EC50 is 1.230-1.519 ng/mL.



ELISA: 5T4 Antibody (13C3) [NBP3-26274] - The Binding Activity of cynomolgus monkey 5T4 with 5T4 Antibody.

Activity: Measured by its binding ability in a functional ELISA.

Immobilized cynomolgus monkey 5T4 at 2 ug/mL can bind 5T4 Antibody, the EC50 is 1.672-2.181 ng/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26274-100ul

| | |
|--------------------|---|
| G-102-C | Goat anti-Human IgG Secondary Antibody [Unconjugated] |
| NB7446 | Goat anti-Human IgG Fc Secondary Antibody |
| DDXCH01P-100 | Human IgG1 Isotype Control |
| H00007162-P01-10ug | Recombinant Human 5T4 GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26274

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

