

Product Datasheet

GIPR Antibody (12C6) - BSA Free NBP3-26553-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26553

Updated 2/24/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26553



NBP3-26553-100ul

GIPR Antibody (12C6) - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20 to -70C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	12C6
Preservative	0.03% Proclin 300
Isotype	IgG2a
Purity	Affinity purified
Buffer	50% Glycerol, 0.01M PBS, PH 7.4

Product Description	
Description	Novus Biologicals Mouse GIPR Antibody (12C6) - BSA Free (NBP3-26553) is a recombinant monoclonal antibody validated for use in ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	2696
Gene Symbol	GIPR
Species	Mouse, Rat
Immunogen	Recombinant Human GIPR protein [UniProt Q0P543]

Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA




Images

ELISA: GIPR Antibody (12C6) [NBP3-26553] - The Binding Activity of Rat GIPR with GIPR Antibody.

Activity: Measured by its binding ability in a functional ELISA.


Immobilized Rat GIPR at 2 ug/mL can bind GIPR Antibody, the EC50 is 6.946-8.740 ng/mL.



ELISA: GIPR Antibody (12C6) [NBP3-26553] - The Binding Activity of Mouse GIPR with GIPR Antibody.

Activity: Measured by its binding ability in a functional ELISA.

Immobilized Mouse GIPR at 2 ug/mL can bind GIPR Antibody, the EC50 is 8.622-11.36 ng/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26553-100ul

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
DBD00	BDNF [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26553

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

