

Product Datasheet

PKHD1 Antibody (8G12A1) - BSA Free NBP3-27118-0.1ml

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27118

Updated 1/19/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27118



NBP3-27118-0.1ml

PKHD1 Antibody (8G12A1) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	8G12A1
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	445 kDa

Product Description	
Description	Novus Biologicals Mouse PKHD1 Antibody (8G12A1) - BSA Free (NBP3-27118) is a monoclonal antibody validated for use in IHC, ELISA, Flow and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rat
Gene ID	5314
Gene Symbol	PKHD1
Species	Human
Immunogen	Purified recombinant fragment of mouse PKHD1 (AA: 3878-4060) expressed in E. Coli.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Flow Cytometry 1:200 - 1:400, ELISA 1:10000, Immunohistochemistry 1:200 - 1:1000, Immunocytochemistry/ Immunofluorescence 1:200 - 1:1000, Immunohistochemistry-Paraffin 1:200 - 1:1000

Images

Flow Cytometry: PKHD1 Antibody (8G12A1) [NBP3-27118] - Analysis of Hela cells using PKHD1 Rat mAb (green) and negative control (red).



Immunohistochemistry: PKHD1 Antibody (8G12A1) [NBP3-27118] - Analysis of paraffin-embedded esophageal cancer tissues using PKHD1 Rat mAb with DAB staining.



ELISA: PKHD1 Antibody (8G12A1) [NBP3-27118] - Black line: Control Antigen (100 ng); Purple line: Antigen (10 ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western Blot: PKHD1 Antibody (8G12A1) [NBP3-27118] - Analysis using PKHD1 mAb against mouse PKHD1 (AA: 3878-4060) recombinant protein. (Expected MW is 23 kDa)



Immunocytochemistry/Immunofluorescence: PKHD1 Antibody (8G12A1) [NBP3-27118] - Analysis of A431 cells using PKHD1 Rat mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27118-0.1ml

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NB7115	Goat anti-Rat IgG (H+L) Secondary Antibody [HRP]
DDXCR01	Rat IgG1 Isotype Control
BC110-60479PEP	PKHD1 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27118

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

