

Product Datasheet

LIR-6/LILRA1 Antibody (4A1E10A6) - BSA Free NBP3-27141-0.1ml

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27141

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27141



NBP3-27141-0.1ml

LIR-6/LILRA1 Antibody (4A1E10A6) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4A1E10A6
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	53.3 kDa

Product Description	
Description	Novus Biologicals Mouse LIR-6/LILRA1 Antibody (4A1E10A6) - BSA Free (NBP3-27141) is a monoclonal antibody validated for use in WB, ELISA and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	11024
Gene Symbol	LILRA1
Species	Human
Immunogen	Purified recombinant fragment of human LIR-6/LILRA1 (AA: extra 322-461) expressed in E. Coli.

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry
Recommended Dilutions	Western Blot 1:500 - 1:2000, Flow Cytometry 1:200 - 1:400, ELISA 1:10000



Images

Flow Cytometry: LIR-6/LILRA1 Antibody (4A1E10A6) [NBP3-27141] - Analysis of MOLT4 cells using LIR-6/LILRA1 mouse mAb (green) and negative control (red).



ELISA: LIR-6/LILRA1 Antibody (4A1E10A6) [NBP3-27141] - Black line: Control Antigen (100 ng); Purple line: Antigen (10 ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western Blot: LIR-6/LILRA1 Antibody (4A1E10A6) [NBP3-27141] - Analysis using LIR-6/LILRA1 mAb against HEK293 (1) and LIR-6/LILRA1 (AA: extra 322-461)-hlgGfc transfected HEK293 (2) cell lysate.



Western Blot: LIR-6/LILRA1 Antibody (4A1E10A6) [NBP3-27141] - Analysis using LIR-6/LILRA1 mAb against human LIR-6/LILRA1 (AA: extra 322-461) recombinant protein. (Expected MW is 41 kDa)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27141-0.1ml

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27141

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

