

Product Datasheet

SP17 Antibody (3G9F3) - BSA Free NBP3-27160-0.1ml

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27160

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27160



NBP3-27160-0.1ml

SP17 Antibody (3G9F3) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3G9F3
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	17.4 kDa

Product Description	
Description	Novus Biologicals Mouse SP17 Antibody (3G9F3) - BSA Free (NBP3-27160) is a monoclonal antibody validated for use in IHC, WB, ELISA and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	53340
Gene Symbol	SPA17
Species	Human
Immunogen	Purified recombinant fragment of human SP17 (AA: 1-152) expressed in E. Coli.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500 - 1:2000, Flow Cytometry 1:200 - 1:400, ELISA 1:10000, Immunohistochemistry 1:200 - 1:1000, Immunohistochemistry-Paraffin 1:200 - 1:1000

Images

Flow Cytometry: SP17 Antibody (3G9F3) [NBP3-27160] - Analysis of SK-OV-3 cells using SP17 mouse mAb (green) and negative control (red).



Western Blot: SP17 Antibody (3G9F3) [NBP3-27160] - Analysis using SP17 mAb against human SP17 (AA: 1-152) recombinant protein. (Expected MW is 20.3 kDa)



Immunohistochemistry: SP17 Antibody (3G9F3) [NBP3-27160] - Analysis of paraffin-embedded rectal cancer tissues using SP17 mouse mAb with DAB staining.



ELISA: SP17 Antibody (3G9F3) [NBP3-27160] - Black line: Control Antigen (100 ng); Purple line: Antigen (10 ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow Cytometry: SP17 Antibody (3G9F3) [NBP3-27160] - Analysis of Hela cells using SP17 mouse mAb (green) and negative control (red).



Immunohistochemistry: SP17 Antibody (3G9F3) [NBP3-27160] - Analysis of paraffin-embedded testis tissues using SP17 mouse mAb with DAB staining.



Western Blot: SP17 Antibody (3G9F3) [NBP3-27160] - Analysis using SP17 mAb against HEK293-6e (1) and SP17 (AA: 1-152)-hlgGfc transfected HEK293-6e (2) cell lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27160-0.1ml

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27160

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

