

Product Datasheet

Proprotein Convertase 9/PCSK9 Antibody [DyLight 650] - (Research Grade tafolecimab Biosimilar) **NBP3-28096C**

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-28096C

Updated 8/14/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-28096C



NBP3-28096C

Proprotein Convertase 9/PCSK9 Antibody [DyLight 650] - (Research Grade tafolecimab Biosimilar)

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2 |
| Conjugate | DyLight 650 |
| Purity | Protein A purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Description | The heavy chain type is hulgG2, and the light chain type is hukappa. It has a predicted MW of 145.5 kDa. Also known as 'tafolecimab'. IgG2SA is human IgG2 with S378A mutation. |
| Host | Human |
| Gene ID | 255738 |
| Gene Symbol | PCSK9 |
| Species | Human |
| Immunogen | PCSK9 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | ELISA, Flow Cytometry, Functional |
| Recommended Dilutions | Flow Cytometry, ELISA, Functional |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-28096C

| | |
|-------------|--|
| NBP3-06870C | Human IgG2 Isotype Control (MAB1641) [DyLight 650] |
| 236-EG-200 | EGF [Unconjugated] |
| MPC900 | Proprotein Convertase 9/PCSK9 [HRP] |
| NB110-60531 | Apolipoprotein E/ApoE Antibody (WUE-4) - BSA Free |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-28096C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

