

# Product Datasheet

## **P-Cadherin Antibody (PF-03732010) [DyLight 650] NBP3-28548C**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-28548C](http://www.novusbio.com/NBP3-28548C)

Updated 8/14/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-28548C](http://www.novusbio.com/reviews/destination/NBP3-28548C)



**NBP3-28548C**

P-Cadherin Antibody (PF-03732010) [DyLight 650]

| <b>Product Information</b> |   |
|----------------------------|---|
| <b>Unit Size</b>           | 0.1 ml  |
| <b>Concentration</b>       | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C in the dark.  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Clone</b>               | PF-03732010   |
| <b>Preservative</b>        | 0.05% Sodium Azide  |
| <b>Isotype</b>             | IgG1  |
| <b>Conjugate</b>           | DyLight 650   |
| <b>Purity</b>              | Protein A purified  |
| <b>Buffer</b>              | 50mM Sodium Borate  |

| <b>Product Description</b> |  |
|----------------------------|--|
| <b>Description</b>         | The heavy chain type is hulgG1, and the light chain type is hulambda. It has a predicted MW of 145.5 kDa.<br><br>Patent #: US20110182884 |
| <b>Host</b>                | Human  |
| <b>Gene ID</b>             | 1001   |
| <b>Gene Symbol</b>         | CDH3   |
| <b>Species</b>             | Human  |
| <b>Immunogen</b>           | CDH3 / P-cadherin  |
| <b>Notes</b>               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.  |

| <b>Product Application Details</b> |  |
|------------------------------------|--|
| <b>Applications</b>                | ELISA, Flow Cytometry, Functional                                      |
| <b>Recommended Dilutions</b>       | Flow Cytometry, ELISA, Functional                                      |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined. |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-28548C**

---

|               |  |
|---------------|--|
| NBP1-85707PEP | P-Cadherin Recombinant Protein Antigen |
| 1129-ER-050   | ErbB2/Her2 [Unconjugated]              |
| 861-PC-100    | P-Cadherin [Unconjugated]              |
| AF748         | E-Cadherin Antibody [Unconjugated]     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-28548C](http://www.novusbio.com/reviews/submit/NBP3-28548C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

