

# Product Datasheet

## **CD11a/CD18 Antibody - Humanized, IgG2SA - Low Endotoxin, Azide and BSA Free NBP3-28609-100ug**

Unit Size: 100 ug

Store at -20C in powder form. Store at -80C once reconstituted.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-28609](http://www.novusbio.com/NBP3-28609)

Updated 12/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-28609](http://www.novusbio.com/reviews/destination/NBP3-28609)



**NBP3-28609-100ug**

CD11a/CD18 Antibody - Humanized, IgG2SA - Low Endotoxin, Azide and BSA Free

| Product Information         |                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ug                                                                                                                                |
| Concentration               | LYOPH mg/ml                                                                                                                           |
| Storage                     | Store at -20C in powder form. Store at -80C once reconstituted.                                                                       |
| Clonality                   | Monoclonal                                                                                                                            |
| Preservative                | No Preservative                                                                                                                       |
| Reconstitution Instructions | Reconstitute with sterile, distilled water to a final concentration of 1 mg/ml. Gently shake to solubilize completely. Do not vortex. |
| Isotype                     | IgG2                                                                                                                                  |
| Purity                      | Protein A purified                                                                                                                    |
| Buffer                      | Lyophilized from 25mM histidine, 8% sucrose, 0.01% Tween80 (pH6.2)                                                                    |

| Product Description |                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | <p>The heavy chain type is hulgG2, and the light chain type is hukappa. It has a predicted MW of 145.5 kDa.</p> <p>This product is shipped at ambient temperature. Upon receipt, store immediately at -20C or lower for 12 months in a lyophilized state. - 80C for 3 months after reconstitution. Avoid repeated freeze-thaw cycles.</p> <p>Patent #: US20160324981</p> |
| Host                | Human                                                                                                                                                                                                                                                                                                                                                                    |
| Gene ID             | 3683                                                                                                                                                                                                                                                                                                                                                                     |
| Gene Symbol         | ITGAL                                                                                                                                                                                                                                                                                                                                                                    |
| Species             | Human                                                                                                                                                                                                                                                                                                                                                                    |
| Immunogen           | Integrin a11 / ITAG11                                                                                                                                                                                                                                                                                                                                                    |
| Endotoxin Note      | < 0.001EU/ug,determined by LAL method.                                                                                                                                                                                                                                                                                                                                   |

| Product Application Details |                                   |
|-----------------------------|-----------------------------------|
| Applications                | ELISA, Flow Cytometry, Functional |
| Recommended Dilutions       | Flow Cytometry, ELISA, Functional |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-28609-100ug**

---

|              |                                                       |
|--------------|-------------------------------------------------------|
| G-102-C      | Goat anti-Human IgG Secondary Antibody [Unconjugated] |
| NB7446       | Goat anti-Human IgG Fc Secondary Antibody             |
| DDXCH02P-100 | Human IgG2 Isotype Control                            |
| 210-TA-005   | TNF-alpha [Unconjugated]                              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-28609](http://www.novusbio.com/reviews/submit/NBP3-28609)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

