

Product Datasheet

Sphingosine Kinase 1/SPHK1 Antibody - BSA Free NBP3-29344

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-29344

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-29344



NBP3-29344

Sphingosine Kinase 1/SPHK1 Antibody - BSA Free

Product Information

Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	Tris-citrate/phosphate buffer, pH 7 to 8

Product Description


Description	Novus Biologicals Rabbit Sphingosine Kinase 1/SPHK1 Antibody - BSA Free (NBP3-29344) is a polyclonal antibody validated for use in IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8877
Gene Symbol	SPHK1
Species	Human
Immunogen	between 1 and 50
Assay Type	[Assay Type]
Sample Volume	[Sample Volume]

Product Application Details

Applications	Immunoprecipitation
Recommended Dilutions	Immunoprecipitation

Images

Immunoprecipitation: Sphingosine Kinase 1/SPHK1 Antibody NBP3-29344] -Detection of human SPHK1 by western blot of immunoprecipitates. Samples: Whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells. Antibodies: Affinity purified rabbit anti-SPHK1 antibody NBP3-29344 used for IP at 10 µg/mg lysate. SPHK1 was also immunoprecipitated by rabbit anti-SPHK1., which recognizes a downstream epitope. For blotting immunoprecipitated SPHK1, was used at 1 µg/ml. Detection: Chemiluminescence with an exposure time of 30 seconds. Note: The more slowly migrating bands are believed to represent additional forms of SPHK1 owing to recognition by multiple antibodies each of which was directed toward distinctly different epitopes on SPHK1.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-29344

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-29344

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

