

Product Datasheet

STAT2 Antibody - BSA Free

NBP3-29475

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-29475

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-29475



NBP3-29475

STAT2 Antibody - BSA Free

Product Information

| | |
|----------------------|---|
| Unit Size | 100 ug |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | Tris-citrate/phosphate buffer, pH 7 to 8 |

Product Description


| | |
|----------------------|---|
| Description | Novus Biologicals Rabbit STAT2 Antibody - BSA Free (NBP3-29475) is a polyclonal antibody validated for use in IHC and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 6773 |
| Gene Symbol | STAT2 |
| Species | Human |
| Immunogen | between 175 and 225 |
| Assay Type | [Assay Type] |
| Sample Volume | [Sample Volume] |

Product Application Details

| | |
|------------------------------|---|
| Applications | Immunohistochemistry, Immunoprecipitation |
| Recommended Dilutions | Immunohistochemistry, Immunoprecipitation |

Images

Immunoprecipitation: STAT2 Antibody NBP3-29475] -Detection of human STAT2 by western blot of immunoprecipitates. Samples: Whole cell lysate (1 mg for IP, 20% of IP loaded) from Jurkat cells. Antibodies: Affinity purified rabbit anti-STAT2 antibody NBP3-29475 used for IP at 6 µg/mg lysate. STAT2 was also immunoprecipitated by rabbit anti-STAT2 antibodies A303-511A and A303-512A, which recognize downstream epitopes. For blotting immunoprecipitated STAT2, A303-512A was used at 1 µg/ml. Detection: Chemiluminescence with an exposure time of 3 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-29475

| | |
|--------------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| H00006773-Q01-10ug | Recombinant Human STAT2 GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-29475

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

