

Product Datasheet

PPAR gamma/NR1C3 Antibody - BSA Free NBP3-30473

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-30473

Updated 1/7/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-30473



NBP3-30473

PPAR gamma/NR1C3 Antibody - BSA Free

Product Information	
Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	Tris-citrate/phosphate buffer, pH 7 to 8

Product Description	
Description	Novus Biologicals Rabbit PPAR gamma/NR1C3 Antibody - BSA Free (NBP3-30473) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5468
Gene Symbol	PPARG
Species	Human, Mouse
Immunogen	The epitope recognized by NBP3-30473 maps to a region between residue 50 to 100 of human Peroxisome Proliferator-Activated Receptor Gamma using the numbering given in entry NP_056953.2 (GeneID 5468).
Assay Type	[Assay Type]
Sample Volume	[Sample Volume]

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:10,000 - 1:25,000, Immunoprecipitation 2 - 10 µg/mg lysate




Images


Immunoprecipitation:[NBP3-30473]-Detection of mouse PPAR-gamma by western blot of immunoprecipitates. Samples: Whole cell lysate from differentiated mouse NIH 3T3-L1 cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit antibody used for IP at 6 µg per reaction. PPAR-gamma was also immunoprecipitated by rabbit antibody. For blotting immunoprecipitated PPAR-gamma, was used at 0. 1 µg/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Whole cell lysate (50 µg) from HeLa, HEK293T, Jurkat, differentiated mouse NIH 3T3-L1, and mouse NIH3T3-L1 cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-PPAR-gamma antibody used for WB at 0.04 µg/ml. Detection: Chemiluminescence with an exposure time of 75 seconds



Whole cell lysate (50 µg) from U937 cells incubated with (+) or without (-) TPA (200nM, 72 hrs) and lysed using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-PPAR-gamma antibody used for WB at 0.1 µg/ml. Detection: Chemiluminescence with an exposure time of 75 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-30473

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-30473

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

