

Product Datasheet

MyoD Antibody (5.2F) - BSA Free NBP3-30815-100ug

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-30815

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-30815



NBP3-30815-100ug


MyoD Antibody (5.2F) - BSA Free

Product Information	
Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5.2F
Preservative	0.08% Sodium Azide
Isotype	IgG2a
Purity	Protein A or G purified
Buffer	Phosphate buffered saline
Product Description	
Description	Novus Biologicals Mouse MyoD Antibody (5.2F) - BSA Free (NBP3-30815) is a monoclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	4654
Gene Symbol	MYOD1
Species	Human, Mouse, Rat, Porcine, Chicken
Immunogen	Recombinant mouse MyoD protein.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot : 1 ug/ml, ELISA, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation : 2 ug/mg of protein lysate, Immunohistochemistry-Paraffin : 2-4 ug/ml, Immunohistochemistry-Frozen : 2-4 ug/ml
Application Notes	.




Images


Immunocytochemistry/ Immunofluorescence: MyoD Antibody (5.2F) [NBP3-30815] - Figure 4. MyoD immunostaining of nuclei in MyoD transfected cells induced to undergo muscle cell differentiation.




Immunocytochemistry/ Immunofluorescence: MyoD Antibody (5.2F) [NBP3-30815] - Figure 3. MyoD immunostaining of nuclei in MyoD transfected cells induced to undergo muscle cell differentiation.



Immunocytochemistry/ Immunofluorescence: MyoD Antibody (5.2F) [NBP3-30815] - Figure 2. MyoD immunostaining of nuclei in MyoD transfected cells induced to undergo muscle cell differentiation.



Western Blot: MyoD Antibody (5.2F) [NBP3-30815] - Left: Immunohistochemical staining using MyoD antibody on formalin fixed, paraffin embedded human rhabdomyosarcoma. Right: Western blot using 60 ng of purified MyoD protein with MyoD antibody at 2 ug/ml and detected using goat anti-mouse HRP and visualized with Pierce West-Femto substrate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-30815-100ug

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-30815

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

