

# Product Datasheet

## Sphingosine Kinase 2/SPHK2 Antibody - BSA Free NBP3-30886-100ug

Unit Size: 100 ug

Store at -70C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-30886](http://www.novusbio.com/NBP3-30886)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-30886](http://www.novusbio.com/reviews/destination/NBP3-30886)



**NBP3-30886-100ug**

Sphingosine Kinase 2/SPHK2 Antibody - BSA Free


Product Information	
Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -70C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.08% Sodium Azide
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	Phosphate buffered saline

Product Description	
Description	Novus Biologicals Rabbit Sphingosine Kinase 2/SPHK2 Antibody - BSA Free (NBP3-30886) is a polyclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	56848
Gene Symbol	SPHK2
Species	Human
Immunogen	Synthetic peptide derived from the Sphingosine Kinase 2/SPHK2 protein.

Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot : 1-5 ug/ml, ELISA
Application Notes	.

**Images**

Western Blot: Sphingosine Kinase 2/SPHK2 Antibody [NBP3-30886] - Western blot of endogenous sphingosine kinase 2 LF detection in Jurkat cell lysate (4 ug/lane) using NBP3-30886 (1 ug/ml) developed using anti-rabbit HRp (1:75K) and Pierce Super Signal west Femto.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-30886-100ug**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-30886](http://www.novusbio.com/reviews/submit/NBP3-30886)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

