

# Product Datasheet

## Cofilin Antibody (PSH0-38)

### NBP3-32206

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-32206](http://www.novusbio.com/NBP3-32206)

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-32206](http://www.novusbio.com/reviews/destination/NBP3-32206)



**NBP3-32206**

Cofilin Antibody (PSH0-38)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PSH0-38
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH7.4), 0.1% BSA and 40% Glycerol
Target Molecular Weight	19 kDa
Product Description	
Description	Novus Biologicals Rabbit Cofilin Antibody (PSH0-38) (NBP3-32206) is a recombinant monoclonal antibody validated for use in WB and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1072
Gene Symbol	CFL1
Species	Human, Mouse, Rat, Monkey
Immunogen	Synthetic peptide within human Cofilin aa 117-166 / 166. (Uniprot: P23528)
Product Application Details	
Applications	Western Blot, Flow Cytometry
Recommended Dilutions	Western Blot 1:1000, Flow Cytometry 1:200-1:1000



## Images

Western Blot: Cofilin Antibody (PSH0-38) [NBP3-32206] - Western blot analysis of Cofilin on different lysates with Rabbit anti-Cofilin antibody (NBP3-32206) at 1/1,000 dilution.

Lane 1: HeLa cell lysate (30 ug/Lane)  
 Lane 2: HEK-293 cell lysate (30 ug/Lane)  
 Lane 3: MCF7 cell lysate (30 ug/Lane)  
 Lane 4: MDA-MB-468 cell lysate (24 ug/Lane)  
 Lane 5: K-562 cell lysate (16 ug/Lane)  
 Lane 6: SH-SY5Y cell lysate (22 ug/Lane)  
 Lane 7: HUVEC cell lysate (30 ug/Lane)  
 Lane 8: Jurkat cell lysate (30 ug/Lane)  
 Lane 9: COS-1 cell lysate (30 ug/Lane)  
 Lane 10: VERO cell lysate (30 ug/Lane)  
 Lane 11: NIH/3T3 cell lysate (22 ug/Lane)  
 Lane 12: PC-12 cell lysate (30 ug/Lane)  
 Lane 13: Neuro-2a cell lysate (19 ug/Lane)

Predicted band size: 19 kDa  
 Observed band size: 19 kDa



Exposure time: 22 seconds;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDN/TBST for 1 hour at room temperature. The primary antibody (NBP3-32206) at 1/1,000 dilution was used in 5% NFDN/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:200,000 dilution was used for 1 hour at room temperature.

Flow Cytometry: Cofilin Antibody (PSH0-38) [NBP3-32206] - Flow cytometric analysis of HeLa cells labeling Cofilin.

Cells were fixed and permeabilized. Then stained with the primary antibody (NBP3-32206, 1ug/ml) (red) compared with Rabbit IgG Isotype Control (green). After incubation of the primary antibody at +4°C for an hour, the cells were stained with a iFluor™ 488 conjugate-Goat anti-Rabbit IgG Secondary antibody at 1/1,000 dilution for 30 minutes at +4°C. Unlabelled sample was used as a control (cells without incubation with primary antibody; black).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-32206**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-32206](http://www.novusbio.com/reviews/submit/NBP3-32206)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

